

Question	Answer	Related Question 1	Related Question 2	Related Question 3				
What is NerveSpa?	NerveSpa is a clinical-grade system that supports patients with neuropathy, joint pain, and circulation issues. It combines aquatic electrotherapy, wraps, lasers, vibration therapy, supplements, and topicals. Products are manufactured in FDA-registered facilities and provided through licensed clinicians.	How do providers prescribe or imple	What conditions and symptoms can Nen	How can my clinic order NerveSpa or get a demo?				
How do providers prescribe or implement NerveSpa?	NerveSpa does not require a medical prescription. It is typically recommended and distributed by licensed providers such as chiropractors, who guide patients on use and may include supplements or consumables as part of a treatment plan.	What is NerveSpa?	What are NerveSpa's pricing options?	How can my clinic order NerveSpa or get a demo?				
What conditions and symptoms can NerveSpa support?	NerveSpa is used for peripheral neuropathy, diabetic neuropathy, chemotherapy-related neuropathy, arthritis (knee, shoulder, OA/RA), plantar fasciitis, heel spurs, restless legs, poor circulation, numbness, tingling, burning, stiffness, balance issues, cramps, cold feet, and general nerve discomfort.	Does NerveSpa help with peripheral	Does NerveSpa help with diabetic neuro	Can NerveSpa help with chemotherapy-induced neuropathy?				
How can my clinic order NerveSpa or get a demo?	Clinics can open a provider account with NerveSpa to place orders or request demos. Starter kits and bundles are available. Patients cannot purchase hardware directly but may obtain consumables through approved channels such as DoctorStore.com.	How do I set up a provider account?	Are starter kits or bundles available for c	How do clinics contact NerveSpa for support?				
What are NerveSpa's pricing options?	Pricing includes wholesale rates for providers (with options for bundles and starter kits) and retail pricing for patients on consumables only, such as salts, garments, and lead wires. Some devices may qualify for insurance coverage under HCPCS code E0720 depending on the plan.	What is the difference between retail	Does insurance ever cover NerveSpa?	Are bundles available for patients?				
What makes NerveSpa unique?	NerveSpa integrates multiple clinical-grade modalities—aquatic electrotherapy, LED wraps, cold laser, vibration therapy—with supplements and topicals to support nerve and joint comfort. Its water-based stimulation provides even, circumferential coverage that pad-only systems can't match.	Is NerveSpa FDA-registered?	How does the Hand & Foot Neuropathy	How does NerveSpa compare to other nerve therapy systems?				
Is NerveSpa FDA-registered?	NerveSpa products are manufactured in FDA-registered facilities and distributed by PMT, an FDA-registered establishment. FDA registration confirms establishment and product listing; it is not an FDA clearance or approval of specific claims.	Is NerveSpa covered by insurance?	Where can FDA listings be verified?	What certifications does PMT hold?				
Is NerveSpa covered by insurance?	Coverage varies by plan. Some providers bill eligible treatments under neuropathy programs; others use cash-pay, HSA/FSA, or memberships. Patients should check benefits with their insurer; clinic billing policies vary.	What is HCPCS code E0720?	Does insurance cover consumables like	How do providers bill insurance for NerveSpa?				
What warranty and return policy does NerveSpa offer? Does NerveSpa replace medical care?	Devices include a limited warranty against defects for one year. Returns are accepted within 30 days if undamaged; buyers are responsible for shipping and a restocking fee. Contact your seller for RMA instructions.	What is the NerveSpa return policy?	Who should patients contact for technical	How quickly are products shipped?				
What is the NerveSpa Hand & Foot Neuropathy System?	No, NerveSpa does not replace medical diagnosis or treatment. Providers determine suitability and protocols; patients should follow provider guidance and report concerns.	Who decides if NerveSpa is appropri	Can NerveSpa be used with other treatm	What role do patients play in outcomes?				
How does the Hand & Foot Neuropathy System work?	The Hand & Foot Neuropathy System is a touch-screen device that uses aquatic electrotherapy with a proprietary waveform. Patients submerge hands or feet in warm water with salts and tablets, allowing circumferential stimulation that supports nerve activity and circulation.	How does the Hand & Foot Neuropa	How long is a typical session?	How often should the system be used?				
How long is a typical session?	The system delivers comfortable stimulation via a warm foot bath to support circulation and nerve health; select Hands/Feet mode and adjust intensity to a mild, comfortable tingle. 30-minute auto-shutoff. Manual: https://drive.google.com/file/d/1ABWl8aU_hrvUSoRkmb3GoHl-rhbhTXlB/view?usp=sharing The default program protocol for 30 minutes and the device auto-shuts off. Dry skin thoroughly after use. Quick Guide: https://drive.google.com/file/d/1ABWl8aU_hrvUSoRkmb3GoHl-rhbhTXlB/view?usp=sharing	Why water-based therapy?	What consumables are used with the sys	What outcomes are expected from the Hand & Foot Neuropathy System?				
How often should the system be used?	Typical clinic/home protocol is once daily, 4–6 days/week for 60–90 days, with continued use as needed. Full User Guide: https://drive.google.com/file/d/1IEghxBYoDlWv3U4V9CDd7D4eulpXcgs/view?usp=sharing	How often should the system be use	Can sessions be longer than 30 minutes?	What happens after the initial 60–90 day program?				
How should the Hand & Foot Neuropathy System be set up?	Connect lead wires and electrodes, add 1 tbsp Epsom salt + 1 effervescent tablet per side, fill with warm water, insert feet, power on, select Hands/Feet, increase intensity until comfortable. Quick Guide: https://drive.google.com/file/d/1ABWl8aU_hrvUSoRkmb3GoHl-rhbhTXlB/view?usp=sharing Full User Guide: https://drive.google.com/file/d/1IEghxBYoDlWv3U4V9CDd7D4eulpXcgs/view?usp=sharing The Hand & Foot Neuropathy System. Charge fully before first use, recharge when low, and avoid leaving it plugged in continuously. A full charge provides up to 72 hours of use.	What improvements do patients rep	What happens after the initial 60–90 day	Can NerveSpa be used long-term?				
How should the battery be maintained?	The device runs a proprietary program with a monophasic waveform and biphasic rest periods at 7.83Hz. Modes automatically alternate to optimize stimulation—patients only adjust intensity.	What if the Hand & Foot Neuropathy	What if stimulation feels weak?	What consumables are used with the system?				
What modes does the system use?	Submersion in warm water allows full circumferential stimulation of the feet or hands. This ensures even current delivery, improved comfort, and deeper engagement of nerves compared to pad-based devices that treat only small surface areas.	What if the system does not power o	How quickly are products shipped?	Who should patients contact for technical issues?				
Why water-based therapy?	Patients often report reduced tingling, burning, and numbness; improved circulation; less stiffness; and better sleep and mobility over time. Outcomes vary by individual and compliance.	Why water-based therapy?	How does this differ from pad-based dev	What outcomes are expected from the Hand & Foot Neuropathy System?				
What improvements do patients report?	Epsom salt and effervescent tablets per side per session; replace as needed. Quick Guide: https://drive.google.com/file/d/1ABWl8aU_hrvUSoRkmb3GoHl-rhbhTXlB/view?usp=sharing NerveSpa Wraps are conductive garment systems designed for localized stimulation of the legs, arms, or joints. They connect to the NerveSpa controller to deliver targeted electrotherapy.	How does this differ from pad-based	What outcomes are expected from the H	How does NerveSpa compare to other nerve therapy systems?				
What consumables are used with the system?	Patients wear the Wraps over the targeted area, moisten them slightly, and connect them to the NerveSpa unit. Providers guide placement to ensure safe and effective stimulation.	Does NerveSpa help with tingling in	Does NerveSpa help with burning sensa	Does NerveSpa help with numbness in the feet or hands?				
What are NerveSpa Wraps?	Patients with pacemakers, other implanted electronics, or severe skin conditions should not use Wraps. Providers should screen patients before recommending.	How do patients reorder consumable	Where can patients get NerveSpa consu	How do clinics order or reorder supplies?				
How are Wraps used in therapy?	The NerveBeam is an advanced OA/RA treatment device using targeted stimulation and wrap design to reduce pain, support cartilage health, and improve knee mobility.	How are Wraps used in therapy?	Who should not use Wraps?	What outcomes are expected from Wraps?				
Who should not use Wraps?	The wrap is positioned around the target area (such as knee, shoulder, or arm), secured with straps, and connected to the controller. Sessions typically last 20–30 minutes as guided by providers.	Who should not use Wraps?	What outcomes are expected from Wrap	Can Wraps be used alongside the Hand & Foot Neuropathy System?				
What is the NerveBeam LED Wrap used?	Patients with pacemakers, other implanted electronics, or light sensitivity conditions should not use the LED Wrap. Providers should screen for contraindications before use.	What precautions should patients fo	Can Wraps be used alongside the Hand	What outcomes are expected from Wraps?				
Who should not use the NerveBeam LED Wrap?	The NerveBeam Cold Laser uses low-level red and infrared laser therapy to reduce inflammation, promote circulation, and support nerve recovery. It is designed for targeted treatment of neuropathy and joint pain.	How is the NerveBeam LED Wrap use	What outcomes are expected from the LED	What outcomes are expected from the LED Wrap?				
Who should not use the NerveBeam Cold Laser?	The handheld applicator is placed directly over the treatment area for 15–20 minutes per session. Providers guide recommended frequency and safety precautions.	Who should not use the Cold Laser?	What outcomes are expected with Cold	Can Cold Laser be used with other NerveSpa devices?				
How is the Cold Laser used?	Patients with pacemakers, epilepsy, or cancerous lesions should not use the Cold Laser. Providers should screen patients for contraindications before use.	Who should not use the Cold Laser?	What outcomes are expected with Cold	Can Cold Laser be combined with other therapies?				
Who should not use the Cold Laser?	The PowerWrap is a high-powered LED and infrared wrap that delivers deep tissue relief for neuropathy, joint pain, and inflammation. It includes 78 lasers and 192 LEDs for clinical-level therapy.	How is the PowerWrap used?	What outcomes are expected with Power	Can the PowerWrap be used with other NerveSpa devices?				
What is the PowerWrap?	The wrap is secured around the target area and connected to its controller. Sessions typically last 20–30 minutes. Providers guide frequency and treatment plans.	Who should not use the PowerWrap?	What outcomes are expected with Powe	Can the PowerWrap be combined with other therapies?				
How is the PowerWrap used?	Patients with pacemakers, other implanted electronics, or certain uncontrolled medical conditions should not use the PowerWrap. Providers should screen before recommending use.	What precautions should patients fo	What outcomes are expected with Powe	Can the PowerWrap be combined with other therapies?				
Who should not use the PowerWrap?	Patients typically report reduced pain, improved circulation, and better mobility over consistent sessions. Outcomes vary by condition and adherence to the program.	What outcomes are expected with th	What outcomes are expected with the K	What outcomes are expected with the ShoulderPro?				
What outcomes are expected with Cold Laser or PowerWrap therapy?	Patients typically report reduced pain, improved mobility, and better circulation when using Wraps consistently. Outcomes vary depending on condition and adherence.	What outcomes are expected with th	What outcomes are expected with the P	What outcomes are expected with the QuakePlate?				
What outcomes are expected with Wrap therapy?	The QuakePlate is a vibrational therapy platform designed to stimulate circulation, relieve pain, and support recovery through targeted vibration therapy.	How is the QuakePlate used?	Who should not use the QuakePlate?	What outcomes are expected with the QuakePlate?				
What is the NerveSpa QuakePlate?	Patients place their feet or hands on the platform and run sessions for 10–20 minutes. Providers can adjust settings based on therapeutic goals.	Who should not use the QuakePlate?	What outcomes are expected with the Q	Can the QuakePlate be used with other NerveSpa devices?				
How is the QuakePlate used?	Patients with pacemakers, seizure disorders, or uncontrolled cardiovascular conditions should not use the QuakePlate. Providers should screen before prescribing.	What precautions should patients fo	What outcomes are expected with the Q	Can the QuakePlate be combined with other therapies?				
Who should not use the QuakePlate?	Patients often report improved circulation, reduced stiffness, and increased comfort in feet, legs, and hands. Results depend on consistency and condition.	What outcomes are expected with th	What outcomes are expected with the S	What outcomes are expected with the Vibe?				
What outcomes are expected with the QuakePlate?	The Kneepro is an advanced OA/RA treatment device using targeted stimulation and wrap design to reduce pain, support cartilage health, and improve knee mobility.	How is the Kneepro used?	Who should not use the Kneepro?	What outcomes are expected with the Kneepro?				
Who should not use the Kneepro?	Patients with pacemakers, open wounds near the knee, or uncontrolled medical conditions should not use the Kneepro. Providers should evaluate safety before prescribing.	Who should not use the Kneepro?	What outcomes are expected with the K	Can the Kneepro be used with other NerveSpa devices?				
What outcomes are expected with the Kneepro?	Patients report reduced knee pain, improved flexibility, and greater mobility with consistent sessions. Outcomes vary depending on severity and adherence.	What precautions should patients fo	What outcomes are expected with the K	Can the Kneepro be combined with other therapies?				
What is the NerveSpa ShoulderPro?	The ShoulderPro is a targeted OA/RA treatment device that delivers therapeutic stimulation to relieve pain, improve mobility, and support joint function in the shoulder.	What outcomes are expected with th	What outcomes are expected with the V	What outcomes are expected with the Hand & Foot Neuropathy System?				
Who should not use the ShoulderPro?	The ShoulderPro is strapped securely around the shoulder and connected to the controller. Sessions generally last 20–30 minutes as recommended by providers.	How is the ShoulderPro used?	Who should not use the ShoulderPro?	What outcomes are expected with the ShoulderPro?				
What outcomes are expected with the ShoulderPro?	Patients with pacemakers, severe skin conditions, or open wounds near the shoulder should not use the ShoulderPro. Providers should evaluate safety before recommending.	Who should not use the ShoulderPro?	What outcomes are expected with the S	Can the ShoulderPro be combined with other therapies?				
What outcomes are expected with the ShoulderPro?	Patients report reduced pain, increased shoulder mobility, and improved daily function with consistent use. Results vary depending on the individual's condition.	What outcomes are expected with th	What outcomes are expected with the H	What outcomes are expected with Wraps?				
What is the NerveSpa Vibe?	The NerveSpa Vibe is a handheld deep tissue vibrational massage designed to relieve muscle tension, support circulation, and complement neuropathy and joint therapies.	How is the Vibe used?	Who should not use the Vibe?	What outcomes are expected with the Vibe?				
Who should not use the Vibe?	It is best applied directly to muscles or joints using different attachment heads. Sessions typically last 10–15 minutes, guided by provider or patient need.	What precautions should patients fo	What outcomes are expected with the V	Can the Vibe be combined with other therapies?				
What outcomes are expected with the Vibe?	Patients with pacemakers, seizure disorders, or severe musculoskeletal conditions should avoid using the Vibe. Providers should screen before recommending.	What outcomes are expected with th	What outcomes are expected with the W	What outcomes are expected with Topicals?				
What topical products are part of NerveSpa?	NerveSpa offers a range of topicals, including Nerve & Neuropathy Cream, Osteoarthritis & RA Cream, Plantar Fasciitis & Heel Spur Cream, and Nerve Target Roll-On.	What is the Nerve & Neuropathy Cream?	What is the Osteoarthritis & RA Cream?	What is the Plantar Fasciitis & Heel Spur Cream?				
What is the Nerve & Neuropathy Cream?	A doctor-formulated topical designed to relieve nerve pain and improve circulation in hands, feet, legs, and joints. Key ingredients include L-Arginine, Vitamin B6, Menthol, and Aloe.	How should the Nerve & Neuropathy Cream be applied?	What is the Plantar Fasciitis & RA Cream?	What is the Osteoarthritis & RA Cream?				
How should the Nerve & Neuropathy Cream be applied?	Apply directly to affected areas such as feet, hands, or joints. Massage gently until absorbed. Use as directed by a provider.	What outcomes are expected with th	What is the Plantar Fasciitis & RA Cream?	What is the Osteoarthritis & RA Cream?				
Who should not use the Nerve & Neuropathy Cream?	Patients with allergies to any listed ingredients, or with severe open wounds on the application area, should avoid using the cream. Consult a provider if unsure.	What outcomes are expected with th	What is the Plantar Fasciitis & RA Cream?	What is the Osteoarthritis & RA Cream?				
What outcomes are expected with the Nerve & Neuropathy Cream?	Patients typically experience reduced nerve pain, improved circulation, and relief from stiffness or inflammation in the applied areas. Results may vary based on condition.	What outcomes are expected with th	What is the Plantar Fasciitis & RA Cream?	What is the Osteoarthritis & RA Cream?				
What outcomes are expected with the Nerve & Neuropathy Cream?	Patients with allergies to any listed ingredients, or with severe skin conditions/open wounds on the application area, should avoid use. Provider consultation recommended.	How should the OA/RA Cream be applied?	What is the Plantar Fasciitis & RA Cream?	What is the Osteoarthritis & RA Cream?				
How should the OA/RA Cream be applied?	Patients often report decreased joint stiffness, reduced pain, and improved mobility. Outcomes depend on consistent use and severity of arthritis.	Who should not use the OA/RA Cream?	What is the Plantar Fasciitis & RA Cream?	What is the Osteoarthritis & RA Cream?				
Who should not use the OA/RA Cream?	A doctor-formulated cream targeting pain and inflammation caused by plantar fasciitis and heel spurs. Supports circulation and provides relief for feet and heels.	What outcomes are expected with th	What is the Plantar Fasciitis & RA Cream?	What is the Osteoarthritis & RA Cream?				
What outcomes are expected with the OA/RA Cream?	Apply to heels or affected foot areas. Massage thoroughly until absorbed. Use consistently as directed by a provider.	How should the OA/RA Cream be applied?	What is the Plantar Fasciitis & RA Cream?	What is the Osteoarthritis & RA Cream?				
What is the Plantar Fasciitis & Heel Spur Cream?	A topical roll-on formula designed to relieve neuropathy symptoms, arthritis, back pain, and sore muscles. Key ingredients include Arnica, Menthol, Camphor, Aloe, and Magnesium.	Who should not use the OA/RA Cream?	What is the Plantar Fasciitis & RA Cream?	What is the Osteoarthritis & RA Cream?				
How should Plantar Fasciitis Cream be applied?	Apply directly to painful areas such as joints, muscles, or feet. Roll on evenly and allow to absorb. Use as directed by a provider.	What outcomes are expected with th	What is the Plantar Fasciitis & RA Cream?	What is the Osteoarthritis & RA Cream?				
What is the Nerve Target Roll-On Pain Relief?	Patients with allergies to Arnica, Menthol, Camphor, or Aloe should avoid use. Avoid applying to broken skin or severe wounds.	How should the OA/RA Cream be applied?	What is the Plantar Fasciitis & RA Cream?	What is the Osteoarthritis & RA Cream?				
How should the Roll-On Pain Relief be applied?	Patients often report reduced muscle tension, relief from inflammation, and cooling/soothing sensations. Relief may be temporary but supports overall therapy.	Who should not use the OA/RA Cream?	What is the Plantar Fasciitis & RA Cream?	What is the Osteoarthritis & RA Cream?				
Who should not use the Roll-On Pain Relief?	Patients typically experience localized relief of pain, inflammation reduction, and improved comfort when used consistently. Results vary by product and condition.	What outcomes are expected with th	What is the Plantar Fasciitis & RA Cream?	What is the Osteoarthritis & RA Cream?				
What outcomes are expected with the Roll-On Pain Relief?	Most patients can safely use NerveSpa topicals, but those with ingredient allergies or severe skin conditions should avoid them. Provider consultation is recommended.	How should the OA/RA Cream be applied?	What is the Plantar Fasciitis & RA Cream?	What is the Osteoarthritis & RA Cream?				
What outcomes are expected from the QuakePlate?	Apply directly to affected areas, massage gently, and allow to absorb. Frequency of use depends on product instructions or provider guidance.	Who should not use the OA/RA Cream?	What is the Plantar Fasciitis & RA Cream?	What is the Osteoarthritis & RA Cream?				
Do topicals replace medical care?	No, NerveSpa creams and roll-ons do not replace medical care. They are supportive tools designed to complement therapies and provider guidance.	What outcomes are expected with th	What is the Plantar Fasciitis & RA Cream?	What is the Osteoarthritis & RA Cream?				
How do topicals integrate with the full NerveSpa system?	Topicals can be used alongside devices and supplements to provide localized relief while supporting systemic improvements.	How should the OA/RA Cream be applied?	What is the Plantar Fasciitis & RA Cream?	What is the Osteoarthritis & RA Cream?				
How do providers determine which topical to recommend?	Providers base recommendations on patient symptoms, conditions, and treatment goals. For example, OA/RA Cream for arthritis or Nerve & Neuropathy Cream for neuropathic pain.	Who should not use the OA/RA Cream?	What is the Plantar Fasciitis & RA Cream?	What is the Osteoarthritis & RA Cream?				
Are demo resources available for creams and roll-ons?	Clinics may receive sample sizes or demo kits to show patients how products are applied and how they complement other NerveSpa therapies.	What outcomes are expected with th	What is the Plantar Fasciitis & RA Cream?	What is the Osteoarthritis & RA Cream?				
Can topicals be used in clinics during therapy sessions?	Yes, topicals are frequently used in conjunction with device-based treatments to maximize comfort and effectiveness.	How should the OA/RA Cream be applied?	What is the Plantar Fasciitis & RA Cream?	What is the Osteoarthritis & RA Cream?				
What support resources are available for providers?	NerveSpa offers provider training, demo videos, and access to clinical support staff. Providers can also request direct assistance for protocols and usage.	Who should not use the OA/RA Cream?	What is the Plantar Fasciitis & RA Cream?	What is the Osteoarthritis & RA Cream?				
What resources are available for patients?	Patients have access to educational videos, FAQs, and support from NerveSpa's customer service team. Providers may also supply resources directly.	What outcomes are expected with th	What is the Plantar Fasciitis & RA Cream?	What is the Osteoarthritis & RA Cream?				
How do clinics contact NerveSpa support?	Providers can reach out for support in the following ways: Chatbot – Use the 'Contact Support' option directly in the chatbot. Email – support@nervespa.com Phone – 1-800-239-7880 (available during business hours). Provider Portal – coming soon.	How should the OA/RA Cream be applied?	What is the Plantar Fasciitis & RA Cream?	What is the Osteoarthritis & RA Cream?				
How do patients contact NerveSpa?	Patients can use the chatbot Contact Support option, call customer service, or email support@nervespa.com for assistance with products, consumables, or supplements.	Who should not use the OA/RA Cream?	What is the Plantar Fasciitis & RA Cream?	What is the Osteoarthritis & RA Cream?				
How do clinics order or reorder supplies?	Providers can place orders directly through their NerveSpa account or by contacting their account representative. Reordering consumables and supplements is streamlined online.	What outcomes are expected with th	What is the Plantar Fasciitis & RA Cream?	What is the Osteoarthritis & RA Cream?				
How do patients receive consumables or supplements?	Patients typically receive consumables (salts, firms, garments, supplements) through their provider, who may also order select consumables directly from DoctorStore.com.	How should the OA/RA Cream be applied?	What is the Plantar Fasciitis & RA Cream?	What is the Osteoarthritis & RA Cream?				
What outcomes are expected from the Hand & Foot Neuropathy System?	Many users report improved circulation, reduced tingling, and greater comfort within weeks. Most consistent outcomes appear over 60–90 days of daily use (4–6×/week) with proper session setup and consumables. See Full User Guide: https://drive.google.com/file/d/1IEghxBYoDlWv3U4V9CDd7D4eulpXcgs/view?usp=sharing	Who should not use the OA/RA Cream?	What is the Plantar Fasciitis & RA Cream?	What is the Osteoarthritis & RA Cream?				
What outcomes are expected from the LED Wrap?	Patients often experience relief from pain, reduced inflammation, and improved tissue recovery when used consistently. https://www.youtube.com/watch?v=GlxoIDKnnw	What outcomes are expected with th	What is the Plantar Fasciitis & RA Cream?	What is the Osteoarthritis & RA Cream?				
What outcomes are expected from the Cold Laser?	Patients often report reduced nerve pain, improved sensation, and faster healing in targeted areas. Clinical testimonials support positive results. https://www.youtube.com/shorts/QRF_zfGeMUl	How should the OA/RA Cream be applied?	What is the Plantar Fasciitis & RA Cream?	What is the Osteoarthritis & RA Cream?				
What outcomes are expected from the PowerWrap?	Patients report targeted relief of neuropathy and joint pain, reduction in inflammation, and improved mobility. Testimonials highlight faster recovery when combined with other therapies. https://www.youtube.com/watch?v=GpyYnLrEQE	Who should not use the OA/RA Cream?	What is the Plantar Fasciitis & RA Cream?	What is the Osteoarthritis & RA Cream?				
What outcomes are expected from the QuakePlate?	Patients typically report decreased shoulder pain, improved range of motion, and relief from arthritis-related stiffness when used consistently.	What outcomes are expected with th	What is the Plantar Fasciitis & RA Cream?	What is the Osteoarthritis & RA Cream?				
What outcomes are expected from the Kneepro?		How should the OA/RA Cream be applied?	What is the Plantar Fasciitis & RA Cream?	What is the Osteoarthritis & RA Cream?				
What outcomes are expected from the ShoulderPro?		Who should not use the OA/RA Cream?	What is the Plantar Fasciitis & RA Cream?	What is the Osteoarthritis & RA Cream?				

Question	Answer
What outcomes are expected from the Vib?	Patients often feel muscle relaxation, reduced stiffness, and improved circulation in targeted areas like legs, back, or arms.
What outcomes are expected from supplements?	Patients typically report improved nerve function, reduced discomfort, enhanced joint support, and improved overall wellness when used as part of a multimodal program.
Can supplements be combined with NerveSpa devices?	Yes, supplements are designed to complement device therapies, creating a multimodal approach for better outcomes.
Are supplements safe for everyone?	Most NerveSpa supplements are safe for adults but may not be suitable for children, pregnant or nursing women, or those with certain medical conditions. Patients should consult their provider.
What ingredients are in N1 Nerve+?	N1 Nerve+ features Korean Red Ginseng in oral dissolvable film format, formulated to improve circulation, boost immunity, and support nervous system health.
How does N1 Nerve+ support nerve health?	Korean Red Ginseng in N1 Nerve+ provides antioxidant support, helps reduce fatigue, and may enhance memory and circulation, all important for nerve repair and function.
What makes the N1 Nerve+ oral dissolvable film unique?	The ODF format avoids added sugars, is easy to take on-the-go, and provides efficient absorption for patients who dislike pills, gummies, or powders.
What ingredients are in N1 Skinny Blend?	N1 Skinny Blend includes Chromium Picolinate, Silver Fir Extract, and Gymnema, supporting weight management, glucose control, and metabolism.
How does N1 Skinny Blend support weight management and metabolism?	Chromium Picolinate helps control sugar intake and cravings, Silver Fir Extract supports oxidative stress defense, and Gymnema promotes glucose utilization and regeneration of islet cells.
What ingredients are in N1 Gut Support?	N1 Gut Support includes Lactobacillus, Probiotics, Saururus Chiniensis Extract, and DL-Malic Acid, supporting digestion, nutrient absorption, and colon health.
How does N1 Gut Support improve digestion and nutrient absorption?	Lactobacillus aids in breaking down food and improving absorption, Propolis provides anti-inflammatory benefits, Saururus Chiniensis supports liver function, and DL-Malic Acid stimulates digestive enzymes.
What makes N1 Gut Support unique?	The ODF format, while supporting probiotic benefits, also supports both digestion and immunity in a convenient daily strip.
What ingredients are in the Super Flex Joint Formula?	The Super Flex Joint Formula contains Glucosamine Sulfate, Chondroitin Sulfate, MSM, Boswellia, Hyaluronic Acid, Turmeric, Ginger, Black Pepper, and White Willow Bark.
How does the Super Flex Joint Formula support joint health?	Glucosamine and Chondroitin stimulate cartilage production, MSM reduces inflammation, and Turmeric and Boswellia add anti-inflammatory and antioxidant support for improved mobility.
What makes the Super Flex Joint Formula unique?	Delivered as a drink powder, it allows higher dosages and encourages hydration. The multi-ingredient formula is designed for patients with osteoarthritis or joint pain needing comprehensive support.
What ingredients are in the Blood Flow Super Formula?	The Blood Flow Super Formula contains L-Arginine, L-Citrulline, vitamins, and minerals designed to increase vasodilation, reduce inflammation, and support nerve and vascular health.
How does the Blood Flow Super Formula support healthy blood circulation?	L-Arginine and L-Citrulline boost nitric oxide production for vasodilation, improving blood flow and circulation while supporting nerve repair and reducing inflammation.
What makes the Blood Flow Super Formula unique?	Delivered as a drink powder, it encourages hydration while combining amino acids with high % daily vitamins. Designed to complement other NerveSpa therapies for circulation and nerve support.
What ingredients are in the Nerve Rebuilder supplement?	The Nerve Rebuilder includes amino acids such as L-Arginine, L-Histidine, L-Glutamine, L-Serine, L-Lysine, Acetyl L-Carnitine, and botanicals like Griffonia Extract, formulated to support nerve repair.
How does the Nerve Rebuilder support nerve health and regeneration?	The blend of amino acids and nutrients helps restore neurotransmitter balance, support nerve regeneration, reduce inflammation, and improve sensation in neuropathy patients.
What makes the Nerve Rebuilder unique?	It is based on clinically researched formulations for neuropathy, with a multi-ingredient blend targeting nerve pain and regeneration. Manufactured in a certified, FDA-registered facility.
What ingredients are in the ImmunityGut Super Formula?	This blend includes Vitamin D, Zinc, Beta Glucan, antioxidants, minerals, and adaptogens, designed for gut restoration and immune system support.
How does the ImmunityGut Super Formula support gut and immune health?	Its three-phase approach—Clear, Seal, Rebuild—strengthens gut lining, reduces inflammation, and fortifies immunity with targeted nutrients.
What makes the ImmunityGut Super Formula unique?	Provides a comprehensive, all-in-one powder addressing gut integrity, detoxification, and stress relief with adaptogens and minerals, unlike single-focus supplements.
How are supplements usually prescribed?	Providers typically recommend supplements as daily support, often alongside device therapy, based on patient needs.
Why is compliance important for supplement outcomes?	Consistent daily use of supplements allows the body to build therapeutic levels of nutrients, maximizing benefits and long-term results.
What ingredients are in the Nerve & Neuropathy Cream?	The cream includes L-Arginine, Vitamin B6, Menthol, and Aloe. These ingredients improve circulation, reduce pain, and soothe inflammation in feet, hands, and legs.
How does the Nerve & Neuropathy Cream provide relief?	L-Arginine boosts nitric oxide for circulation, Menthol provides cooling pain relief, and Vitamin B6 supports nerve health. The cream absorbs quickly for targeted action.
What ingredients are in the Osteoarthritis & RA Cream?	The OA/RA Cream contains Glucosamine, Arnica, Arginine, Menthol, and Aloe, designed to reduce inflammation, rebuild cartilage, and relieve joint pain.
How does the Osteoarthritis & RA Cream help?	Glucosamine supports cartilage rebuild, Arnica and Menthol reduce inflammation and pain, and Aloe soothes tissues. The cream absorbs quickly for targeted joint relief.
What is the Nerve Target Roll-On?	Nerve Target Roll-On is a topical pain reliever with Arnica, Menthol, Camphor, and Aloe, designed for neuropathy, arthritis, and muscle pain.
How should the Nerve Target Roll-On be applied?	Apply directly to affected areas such as knees, back, or feet. The roll-on format ensures easy application without mess.
What outcomes are expected from creams and roll-ons?	Patients report reduced localized pain, improved comfort in daily activities, and better circulation when topicals are applied consistently as directed.
Are topicals safe for everyone?	Generally safe, but patients with skin sensitivities or allergies should check ingredients or consult a provider before use.
Can topicals be used with other NerveSpa therapies?	Yes. Creams and roll-ons can be paired with devices and supplements for a multi-modal approach, enhancing overall outcomes.
Do topicals replace oral supplements or devices?	No. Topicals provide localized relief, while devices and supplements address underlying nerve, joint, or circulation issues. They are intended to work together.
Why is compliance important for topical outcomes?	Regular application ensures active ingredients build consistent therapeutic effects, supporting pain relief and circulation benefits.
What support resources exist for topical use?	Clinics may provide application demos, and NerveSpa offers instructions, videos, and customer support for correct use of creams and roll-ons.
How do patients reorder creams and roll-ons?	Patients can reorder through their provider or directly from approved online channels such as doctorstore.com.
Do topicals require FDA clearance?	NerveSpa creams and roll-ons are made in FDA-registered and cGMP facilities, but they are not FDA-cleared medical devices.
Do providers need to prescribe topicals?	Topicals do not require a prescription. Providers typically recommend them as part of a patient's care plan.
Are topicals covered by warranty or return policies?	Yes. NerveSpa offers satisfaction guarantees and return policies for creams and roll-ons. Clinics can confirm details with their representative.
Who is PMT in relation to NerveSpa?	PMT (Pain Management Technologies) is the parent company of NerveSpa, with over 30 years of experience in developing clinical-grade electrotherapy devices.
What certifications does PMT hold?	PMT (Pain Management Technologies) is an FDA-registered establishment. NerveSpa products are produced in FDA-registered facilities, and supplements and topicals are manufactured in FDA-registered and cGMP facilities.
Where are NerveSpa products manufactured?	NerveSpa products are sourced globally. Devices are produced in FDA-registered facilities, while supplements and topicals are made in FDA-registered and cGMP facilities.
Where can FDA listings be verified?	PMT and NerveSpa's FDA establishment and product registrations can be verified in the FDA registration database or via the official spreadsheet: https://bit.ly/41mWufw
What warranty does NerveSpa offer?	NerveSpa devices include a standard limited warranty. Terms vary by product and can be confirmed with a NerveSpa representative.
What is the NerveSpa return policy?	NerveSpa offers a return policy for devices and consumables if patients are not satisfied. Specific terms should be confirmed at purchase or with the provider.
Does NerveSpa replace medical treatment?	No. NerveSpa is a complementary therapy. It does not replace medical care, diagnosis, or treatment from a licensed healthcare professional.
Who determines how NerveSpa should be used?	NerveSpa protocols are typically recommended by providers such as chiropractors or physical therapists. Patients use the devices at home under provider guidance.
Do NerveSpa products require direct provider supervision?	No. Products are designed for home use. Providers prescribe or recommend protocols, but patients can use devices independently.
What should patients do if discomfort or side effects occur?	Stop use immediately and consult their provider. Side effects are rare but any discomfort should be reported to a healthcare professional.
What if the Hand & Foot Neuropathy System does not power on?	Check that the charger is connected, the device is fully charged, and all cables are secure. If issues persist, contact NerveSpa support.
What if stimulation feels weak?	Ensure salt and effervescent tablets are added to the water, electrodes are connected properly, and feet or hands are fully submerged. If still weak, consult support.
What if the LED Wrap does not light?	Confirm that the charger is connected and the wrap is powered on. Check replacement parts if needed. If unresolved, contact NerveSpa support.
How does NerveSpa ensure safety?	All NerveSpa products are manufactured in FDA-registered facilities. Supplements and topicals are produced in FDA-registered and cGMP facilities. Devices are drug-free, non-invasive, and designed for home or clinic use.
What if Wraps feel uncomfortable?	Ensure wraps are secured snugly but not too tight. Adjust position for comfort. Discontinue use if irritation occurs and consult support if needed.
How quickly are products shipped?	Most orders ship within 1-3 business days. Shipping times may vary depending on product availability and destination.
How are returns processed?	Returns are handled through NerveSpa or the clinic where the purchase was made. Contact support for instructions and eligibility.
How do providers bill insurance for NerveSpa?	Insurance coverage is uncommon. Providers should consult payer guidelines for possible reimbursement of consumables or related therapies.
How do clinics contact NerveSpa for support?	Clinics can contact NerveSpa through the official support portal, email, or phone for product or billing questions.
What role do patients play in outcomes?	Patients are key to success. Consistency with therapy, adherence to provider guidance, and healthy lifestyle choices all enhance results.
What makes N1 Nerve+ oral dissolvable film unique?	N1 Nerve+ uses an oral dissolvable film (ODF) delivery system—slim, easy to take on-the-go, and ideal for patients who avoid pills or powders.
Do providers prescribe supplements or can patients self-select?	Supplements are generally recommended by providers for best results, though patients can purchase them directly when appropriate.
What is the Osteoarthritis & Rheumatoid Arthritis (OA/RA) Cream?	The OA/RA Cream combines Glucosamine, Arnica, Arginine, Menthol, and Aloe to relieve joint pain, reduce inflammation, and support cartilage and tissue repair.
How should the OA/RA Cream be used?	Apply directly to painful joints (hands, knees, feet). It absorbs quickly and can be used to support daily joint health alongside therapy.
Where can patients get NerveSpa creams and roll-ons?	Patients can obtain them through their provider or directly from NerveSpa. Consumables may also be available via doctorstore.com.
What support resources are available for topicals?	Support includes training guides, application videos, and direct access to NerveSpa's clinical support team for questions or troubleshooting.
How do patients say about creams and roll-ons?	Patients report fast relief, improved circulation, reduced inflammation, and ease of use with the roll-on format or fast-absorbing creams.
What do patients or clinics order topicals?	Orders can be placed through NerveSpa, authorized providers, or doctorstore.com for consumables. Clinics may set up wholesale accounts for bulk supply.
Are there minimum order requirements?	Yes. NerveSpa offers training guides and educational resources to providers, understand indications, application, and integration with other therapies.
Can topicals be used long-term?	Yes. They are safe for ongoing use when applied as directed, making them a useful part of chronic neuropathy or arthritis care plans.
Does insurance cover consumables like salts or garments?	Consumables (e.g., salts, garments, lead wires) are not typically covered by insurance. Patients usually purchase them directly via providers or doctorstore.com.
What is HCPCS code E0720?	HCPCS code E0720 refers to transcutaneous electrical nerve stimulation (TENS). Some providers may use it when billing electrotherapy treatments, but coverage varies by insurer.
How do clinics order NerveSpa?	Clinics can order by contacting NerveSpa directly to set up a provider account. A wholesale portal is planned but not yet live.
How do clinics bill insurance for NerveSpa?	Clinics may use HCPCS code E0720 for electrotherapy when appropriate. Coverage depends on the patient's insurer and documentation.
How do I set up a provider account?	Providers should contact NerveSpa directly to establish an account. This grants access to wholesale pricing, training, and ordering. A dedicated online portal will be added in the future.
How do providers bill insurance for NerveSpa?	Pricing for patients depends on the products recommended by their provider. Devices, consumables, and supplements can be ordered via providers or doctorstore.com.
What products are included in a starter kit?	Provider accounts may include minimum purchase quantities for wholesale pricing. Requirements vary by product type and order size.
Are bundles available for patients?	Starter kits typically include the Hand & Foot Neuropathy System, consumables (salt, tablets, lead wires), and educational materials. Options vary by clinic needs.
Who should patients contact for technical issues?	Yes. NerveSpa offers patient bundles that may include devices, supplements, and topicals to support therapy at home. Availability may vary by clinic.
Are starter kits or bundles for clinics available?	Patients should contact NerveSpa support through the website contact form or by phone for troubleshooting, warranty claims, or replacement requests.
What outcomes are expected from bundles?	Yes. NerveSpa offers starter kits and clinic bundles designed for easy onboarding. These include core devices, consumables, supplements, and educational materials.
What training do providers receive?	Bundles provide patients and clinics with a complete therapy package. Combining devices, supplements, and topicals helps improve consistency and outcomes.
How can patients learn proper use?	Providers receive training resources including setup guides, demo videos, and optional live onboarding calls to ensure correct use.
What outcomes are expected from Wraps?	Patients can watch demo videos, read included manuals, or request live support from their provider or NerveSpa.
What outcomes are expected from conductive garments?	Wraps help deliver consistent stimulation to targeted areas, improving comfort, reducing pain, and supporting recovery.
Is NerveSpa safe for everyone?	Conductive garments (socks, gloves) help deliver even stimulation across the entire hand or foot, supporting circulation, reducing pain, and calming nerve discomfort.
What support resources are available?	NerveSpa devices and products are designed to be safe when used as directed, but contraindications apply (e.g., pacemakers, open wounds). Providers should guide usage.
Can NerveSpa be used with other treatments?	NerveSpa provides provider guides, demo videos, manuals, and direct support to help with onboarding and patient use.
	Yes. NerveSpa can be safely combined with many therapies such as chiropractic care, PT, or red light therapy. Providers determine the best plan.

Related Question 1	Related Question 2	Related Question 3
What outcomes are expected with to What outcomes are expected with suppl	What outcomes are expected with suppl	What outcomes are expected with the Hand & Foot Neuropathy System?
What ingredients are in N1 Nerve+?	What ingredients are in N1 Skinny Blend	What ingredients are in N1 Gut Support?
What outcomes are expected from s	Are supplements safe for everyone?	Why is compliance important for supplement outcomes?
What outcomes are expected from s	How are supplements usually prescribed	Why is compliance important for supplement outcomes?
How does N1 Nerve+ support nerve	What makes the N1 Nerve+ oral dissolve	What outcomes are expected from supplements?
What makes the N1 Nerve+ oral dissolvable film unique?	What makes the N1 Nerve+ oral diss	What outcomes are expected from supplements?
How does N1 Nerve+ support nerve	What ingredients are in N1 Skinny Blend	What ingredients are in N1 Gut Support?
How does N1 Skinny Blend support	What outcomes are expected from suppl	What makes the N1 Skinny Blend unique?
What makes the N1 Skinny Blend	What ingredients are in N1 Gut Support?	Are supplements safe for everyone?
How does N1 Gut Support improve	What outcomes are expected from suppl	Why is compliance important for supplement outcomes?
What makes N1 Gut Support unique	What outcomes are expected from suppl	Are supplements safe for everyone?
What makes the N1 Gut Support unique?	What makes the N1 Gut Support unique?	What outcomes are expected from supplements?
What makes the Super Flex Joint Formu	What makes the Super Flex Joint Formu	What outcomes are expected from supplements?
What makes the Super Flex Joint Formu	What outcomes are expected from suppl	Are supplements safe for everyone?
What ingredients are in the Blood Fk	How does the Blood Flow Super Formu	What outcomes are expected from supplements?
How does the Blood Flow Super For	What makes the Blood Flow Super Form	What outcomes are expected from supplements?
What makes the Blood Flow Super F	What outcomes are expected from suppl	Are supplements safe for everyone?
What ingredients are in the Nerve R	How does the Nerve Rebuilder support	n
What ingredients are in the Nerve Rebuilder	What outcomes are expected from suppl	What outcomes are expected from supplements?
How does the Nerve Rebuilder supp	What makes the Nerve Rebuilder unique	What outcomes are expected from supplements?
What makes the Nerve Rebuilder un	What outcomes are expected from suppl	Are supplements safe for everyone?
What ingredients are in the ImmunityGut Super Form,	What outcomes are expected from suppl	What outcomes are expected from supplements?
How does the ImmunityGut Super F	What makes the ImmunityGut Super For	What outcomes are expected from supplements?
What makes the ImmunityGut Super	What outcomes are expected from suppl	Are supplements safe for everyone?
What outcomes are expected from s	Can supplements be combined with Nen	Why is compliance important for supplement outcomes?
What outcomes are expected from s	Why is compliance important for supplien	Are supplements safe for everyone?
What outcomes are expected from s	Can supplements be combined with Nen	Are supplements safe for everyone?
How does the Nerve & Neuropathy C	What outcomes are expected from crean	Are topicals safe for everyone?
What ingredients are in the Osteoart	What outcomes are expected from crean	Are topicals safe for everyone?
How does the Osteoarthritis & RA Cr	What outcomes are expected from crean	Are topicals safe for everyone?
What is the Plantar Fasciitis & Heel	What outcomes are expected from crean	Are topicals safe for everyone?
How should the Nerve Target Roll-O	What outcomes are expected from crean	Are topicals safe for everyone?
What outcomes are expected from c	Are topicals safe for everyone?	Can topicals be used with other NerveSpa therapies?
How should NerveSpa creams and r	Can topicals be combined with other Ner	Why is compliance important for topical outcomes?
What outcomes are expected from c	How should NerveSpa creams and roll-o	Can topicals be combined with other NerveSpa therapies?
What outcomes are expected from c	Why is compliance important for topical	c
What outcomes are expected from c	Can topicals be combined with other Ner	Why is compliance important for topical outcomes?
What outcomes are expected from c	Can topicals be combined with other Ner	What support resources exist for topical use?
What outcomes are expected from c	How do patients reorder creams and roll	How do clinics order or reorder supplies?
How do clinics order or reorder supp	What support resources exist for topical	Are topicals covered by warranty or return policies?
Are topicals safe for everyone?	What outcomes are expected from crean	Do topicals replace oral supplements or devices?
How do patients reorder creams and	How do clinics order or reorder supplies	What support resources exist for topical use?
How do patients reorder creams and	How do clinics order or reorder supplies	What support resources exist for topical use?
What certifications does PMT hold?	Where are NerveSpa products manufact	How does NerveSpa ensure safety?
Who is PMT in relation to NerveSpa?	Where are NerveSpa products manufact	Where can FDA listings be verified?
What certifications does PMT hold?	Where can FDA listings be verified?	How does NerveSpa ensure safety?
What certifications does PMT hold?	Is NerveSpa FDA-registered?	Where are NerveSpa products manufactured?
What is the NerveSpa return policy?	Are topicals covered by warranty or retur	Who should patients contact for technical issues?
What warranty does NerveSpa offer?	Are topicals covered by warranty or retur	Who should patients contact for technical issues?
Who determines how NerveSpa sho	What support resources are available?	What role do patients play in outcomes?
Does NerveSpa replace medical tres	What support resources are available?	How do providers monitor patient progress?
Who determines how NerveSpa sho	Does NerveSpa replace medical treatments	What support resources are available?
Is NerveSpa safe for older adults?	Who should not use NerveSpa products	What role do patients play in outcomes?
What if stimulation feels weak?	Who should patients contact for technica	How are returns processed?
What if the Hand & Foot Neuropathy	What if Wraps feel uncomfortable?	Who should patients contact for technical issues?
What if the Cold Laser does not wor	Who should patients contact for technica	Are demo resources available for creams and roll-ons?
What certifications does PMT hold?	Is NerveSpa FDA-registered?	Who determines how NerveSpa should be used?
What if stimulation feels weak?	What if the LED Wrap does not light?	Who should patients contact for technical issues?
How are returns processed?	What warranty does NerveSpa offer?	How do patients receive consumables or supplements?
What warranty does NerveSpa offer?	What is the NerveSpa return policy?	Who should patients contact for technical issues?
Does insurance ever cover NerveSp	Does insurance cover consumables like	How do clinics order NerveSpa?
How do patients contact NerveSpa?	Who should patients contact for technica	What support resources are available?
What lifestyle changes support Nerv	Can NerveSpa be used long-term?	How do providers monitor patient progress?
What ingredients are in N1 Nerve+?	How does N1 Nerve+ support nerve hea	What outcomes are expected from supplements?
Are supplements safe for everyone?	How are supplements usually prescribed	What outcomes are expected from supplements?
How should the OA/RA Cream be us	What outcomes are expected from topici	Are topicals safe for everyone?
What is the Osteoarthritis & RA Crea	What outcomes are expected from topici	Are topicals safe for everyone?
How do patients reorder creams and	Do topicals require FDA clearance?	Do topicals replace medical care?
Are demo resources available for cre	How do patients reorder creams and roll	What do patients say about creams and roll-ons?
What outcomes are expected from N	Are topicals safe for everyone?	How do providers determine which topical to recommend?
Where can patients get creams and	How do patients reorder creams and roll	Are topicals covered by warranty or return policies?
Are demo resources available for cre	What support resources are available for	How do topicals integrate with the full NerveSpa system?
Are topicals safe for everyone?	Do topicals replace medical care?	What outcomes are expected from NerveSpa topicals?
How can patients order consumable	What are NerveSpa's pricing options?	Is NerveSpa covered by insurance?
Is NerveSpa covered by insurance?	Does NerveSpa require a prescription?	How do clinics bill insurance for NerveSpa?
How do I set up a provider account?	Are there minimum order requirements?	Are starter kits or bundles for clinics available?
What is HCPCS code E0720?	Is NerveSpa covered by insurance?	How do clinics or patients contact NerveSpa?
How do clinics order NerveSpa?	Are there minimum order requirements?	Are starter kits or bundles for clinics available?
How do patients order NerveSpa?	Does insurance cover NerveSpa?	How can patients order NerveSpa?
Are bundles available for patients?	How do I set up a provider account?	Are starter kits or bundles for clinics available?
What products are included in a star	How do clinics order NerveSpa?	How do I set up a provider account?
How do clinics contact NerveSpa for	What is the NerveSpa return policy?	What warranty does NerveSpa offer?
What products are included in a star	How do clinics order NerveSpa?	Are there minimum order requirements?
Are bundles available for patients?	Are starter kits or bundles for clinics avai	What resources are available for patients?
What support resources are availabl	How do clinics order NerveSpa?	How can patients learn proper use?
What resources are available for pat	What training do providers receive?	Who should patients contact for technical issues?
What are NerveSpa Wraps?	How are Wraps used in therapy?	Who should not use Wraps?
What are NerveSpa conductive gar	How are conductive garments used?	Who should not use conductive garments?
Who should not use the Hand & Foo	Who should not use Wraps?	Who should not use the Cold Laser?
What training do providers receive?	What resources are available for patients	How can patients learn proper use?
Can multiple NerveSpa modalities b	What lifestyle changes support NerveSp	Who decides if NerveSpa is appropriate?

Question	Answer	Related Question 1	Related Question 2	Related Question 3					
Can NerveSpa be used long-term?	NerveSpa is designed for sustained use. After the initial 60–90 day program, most patients transition to a maintenance schedule recommended by their provider.	What happens after the initial 60–90	What lifestyle changes support NerveSpa?	Who decides if NerveSpa is appropriate?					
Can multiple NerveSpa modalities be used in the same day?	Yes. Patients often combine the Hand & Foot Neuropathy System with wraps, lasers, or topicals for complementary results, under provider guidance.	Can NerveSpa be used with other tr	What outcomes are expected from comb	Who decides if NerveSpa is appropriate?					
What happens after the initial 60–90 day program?	After 60–90 days of consistent use, patients usually shift to a maintenance schedule with reduced frequency, as guided by their provider.	Can NerveSpa be used long-term?	Can NerveSpa be used with other treatm	What outcomes are expected from the Hand & Foot Neuropathy System?					
What lifestyle changes support NerveSpa therapy?	Why by: nutrition, exercise, and limiting alcohol and tobacco can enhance outcomes with NerveSpa therapies.	Can NerveSpa be used long-term?	Who decides if NerveSpa is appropri	What role do patients play in outcomes?					
What research supports NerveSpa therapies?	NerveSpa is built on established electrotherapy, red light therapy, and supplement research. Providers can request clinical references, white papers, and case studies from the NerveSpa team.	Who decides if NerveSpa is appropri	What role do patients play in outcomes?	How does NerveSpa ensure safety?					
Who decides if NerveSpa is appropriate?	Providers evaluate patient history, symptoms, and goals to determine if NerveSpa is suitable, and which modalities or supplements to use.	What research supports NerveSpa t	What role do patients play in outcomes?	What support resources are available?					
What should patients do if discomfort or side effects occur?	Stop use immediately and consult the prescribing provider. If symptoms persist, contact NerveSpa support for guidance.	Is NerveSpa safe for everyone?	Who should not use NerveSpa products?	What support resources are available?					
What support resources are available for patients?	Patients can access demo videos, user manuals, supplement guides, and request live assistance through their provider or NerveSpa support.	How do patients contact NerveSpa?	What resources are available for patients?	How can patients learn proper use?					
Who should not use NerveSpa Wraps?	NerveSpa Wraps should not be used by patients with open wounds, active skin infections, implanted electronic devices (such as pacemakers), or by those who have been specifically advised by their provider not to use electrical stimulation therapies. Wraps are intended for external use only.	If patients want topical options, what	What other NerveSpa devices are availa	Can Wraps be combined with other therapies?					
Who should not use NerveSpa Creams and Roll-Ons?	Creams and Roll-Ons should not be applied to broken skin, open wounds, or areas of active rash or infection. Patients with allergies to listed ingredients (menthol, camphor, arnica, etc.) should avoid use.	Where can patients see ingredients	How do patients reorder creams and roll	Can topicals be used alongside devices?					
Can NerveSpa topicals be combined with device therapies?	Yes. NerveSpa creams and Roll-Ons can be safely combined with device therapies. They are often recommended as part of a multimodal program, providing targeted relief while stimulation supports circulation and nerve activity.	What outcomes are expected from c	What outcomes are expected from devic	What lifestyle changes can support therapy?					
How do patients reorder NerveSpa supplements?	Patients can reorder supplements through their provider's office or approved distributor channels such as doctorstore.com. Providers should guide on which supplement is appropriate for their condition and ensure compliance with the overall program.	How do providers reorder supplie	What outcomes are expected from suppl	What support is available if there are ordering issues?					
How do patients reorder consumables like salts, electrodes, or lead wires?	Consumables such as salts, electrodes, and lead wires are reordered through provider clinics or approved portals like doctorstore.com. Providers may also set up automatic reordering systems for long-term patients.	What consumables are used with the	How quickly are products shipped?	What support is available if there are ordering issues?					
Can NerveSpa supplements be combined with device therapies?	Yes. Supplements are designed to work in tandem with NerveSpa devices by supporting circulation, nutrient supply, and nerve repair at the systemic level while devices target local nerve stimulation. Lifestyle changes such as maintaining balanced nutrition, consistent physical activity, proper hydration, sleep hygiene, and limiting alcohol or tobacco use can significantly enhance NerveSpa therapy outcomes. Stress management and compliance with provider guidance are also important.	What outcomes are expected from s	What outcomes are expected from devic	What lifestyle changes can support therapy?					
What lifestyle changes can support NerveSpa therapy?	Providers can reach out for support in the following ways: Chatbot – Use the "Contact Support" option directly in the chatbot. Email – support@nervespa.com Phone – 1-800-239-7880 (available during business hours) Provider Portal – coming soon.	Can NerveSpa be used long-term?	What role do patients play in outcomes?	What support resources are available?					
Contact phone number email	Phone – 1-800-239-7880 (available during business hours). Email – support@nervespa.com	phone number	email						
Who should NOT use NerveSpa Hand & Footbath?	Do not use with open wounds on treatment area, during pregnancy, or with implanted pacemakers without physician guidance. Check skin frequently. Full User Guide: https://drive.google.com/file/d/1IEghxB0y0alwv3U4v9CDd7D4eulpxcgs/view?usp=sharing	Contact	email						
What are NerveSpa Conductive Garments (gloves/socks)?	Alternative to footbath using conductive gloves/socks with the NerveSpa device. Spray inside/outside until damp, wear comfortably, and adjust intensity to a mild tingle. Guide: https://drive.google.com/file/d/18q1b0FWLzUbsCJmdrwuefh3oEnUSlC/view?usp=sharing Video: https://www.youtube.com/@nervespa	Contact	phone number						
How do I use the conductive gloves or socks?	Charge device, connect lead wires to garment pigtails, spray garments inside and out, put on, start program and adjust intensity as needed. Use 4–6x/week. Guide: https://drive.google.com/file/d/18q1b0FWLzUbsCJmdrwuefh3oEnUSlC/view?usp=sharing Video: https://www.youtube.com/@nervespa	Is NerveSpa safe?	Can seniors use NerveSpa at home?	Is NerveSpa covered by insurance?					
How often should conductive garments be used?	Use once daily, 4–6 days/week. Benefits typically accumulate over 60–90 days. Guide: https://drive.google.com/file/d/18q1b0FWLzUbsCJmdrwuefh3oEnUSlC/view?usp=sharing	How do I use the conductive garmen	How often should the system be used?	How do I clean and store accessories?					
Who should NOT use the conductive garments?	Do not use if under 18, pregnant, have a pacemaker, open wounds, total loss of sensation, or skin irritation; check skin every 15 minutes. Guide: https://drive.google.com/file/d/18q1b0FWLzUbsCJmdrwuefh3oEnUSlC/view?usp=sharing	What are NerveSpa Conductive Can	What does treatment feel like?	Who should NOT use NerveSpa?					
What is NerveBeam (LED/Power Wrap)?	Use once daily, 4–6 days/week. Benefits typically accumulate over 60–90 days. Guide: https://drive.google.com/file/d/18q1b0FWLzUbsCJmdrwuefh3oEnUSlC/view?usp=sharing	What does treatment feel like?	How long is a typical session?	How do I clean and store accessories?					
How do I use the NerveBeam LED/Power Wrap?	Do not use if under 18, pregnant, have a pacemaker, open wounds, total loss of sensation, or skin irritation; check skin every 15 minutes. Guide: https://drive.google.com/file/d/18q1b0FWLzUbsCJmdrwuefh3oEnUSlC/view?usp=sharing	Who should NOT use NerveSpa Har	Is NerveSpa safe?	What does treatment feel like?					
How long is a NerveBeam session?	Flexible red (660nm) and infrared (830/850nm) light therapy wraps for neuropathy, pain, and recovery. LED Wrap Manual: https://drive.google.com/file/d/11FLZ7E108Iyx4SdKSiBHP-NNskbOmC/view?usp=sharing Power Wrap Manual: https://drive.google.com/file/d/1tWDJ7dD9kanGU925nk0Kp-17t09aupF5/view?usp=sharing Video: https://www.youtube.com/@nervespa	How do I use the NerveBeam wrap?	What are the benefits of red vs infrared?	How long is a NerveBeam session?					
Who should NOT use NerveBeam?	Plug in or charge, position wrap over area, secure LOOSELY (allow airflow), press power 2s. Red only for skin; red+IR for deeper tissues. Auto-shutoff ~20 min. LED Manual: https://drive.google.com/file/d/11FLZ7E108Iyx4SdKSiBHP-NNskbOmC/view?usp=sharing Power Wrap Manual: https://drive.google.com/file/d/1tWDJ7dD9kanGU925nk0Kp-17t09aupF5/view?usp=sharing Video: https://www.youtube.com/@nervespa	What is NerveBeam (LED/Power Wr	What does treatment feel like?	Who should NOT use NerveBeam?					
How do I clean and store the NerveBeam?	Avoid use on malignant tumors, open wounds, recent injuries (<72h), and avoid tight wrapping/heat trapping. Not for pregnancy; supervise elderly/children; do not sleep with device. LED Manual: https://drive.google.com/file/d/11FLZ7E108Iyx4SdKSiBHP-NNskbOmC/view?usp=sharing Power Wrap Manual: https://drive.google.com/file/d/1tWDJ7dD9kanGU925nk0Kp-17t09aupF5/view?usp=sharing	How do I use the NerveBeam LED/P	What does treatment feel like?	Who should NOT use NerveBeam?					
What is QuakePlate and what does it do?	A deep-tissue vibration/massage platform for feet/lower extremities to improve mobility, relieve pain, and support circulation. Manual: https://drive.google.com/file/d/1BNgb7AoSEeZTWMBtTh5R7r5xXs1b_n/view?usp=sharing	How do I use the NerveBeam LED/P	Who should NOT use NerveBeam?	How long is a NerveBeam session?					
How do I use QuakePlate safely?	Use indoors on 110V AC, keep connections dry, unplug after use, do not use with wet hands. Monitor comfort; stop with any adverse reaction. Manual: https://drive.google.com/file/d/1BNgb7AoSEeZTWMBtTh5R7r5xXs1b_n/view?usp=sharing	How do I use QuakePlate safely?	Who should NOT use QuakePlate?	How long is a QuakePlate session?					
Who should NOT use QuakePlate?	Avoid if pregnant, with recent injuries, local infections, or circulatory syndromes (e.g., Raynaud's, PVD, sickle cell). Manual: https://drive.google.com/file/d/1BNgb7AoSEeZTWMBtTh5R7r5xXs1b_n/view?usp=sharing	What is QuakePlate and what does i	Who should NOT use QuakePlate?	How do I clean and maintain QuakePlate?					
How do I clean and maintain QuakePlate?	Wipe surfaces, check cords away from heat, inspect for damage; unplug immediately if cords are damaged; contact supplier for service. Manual: https://drive.google.com/file/d/1BNgb7AoSEeZTWMBtTh5R7r5xXs1b_n/view?usp=sharing	Is QuakePlate safe?	How do I use QuakePlate safely?	How long is a QuakePlate session?					
What is NerveSpa Knee Pro?	A knee wrap with Jstlm electrotherapy for OA/RA therapy aimed at regenerative support (cartilage, tendons, ligaments). Quick Start: https://drive.google.com/file/d/17HlOn5mNzZxvJvtuPjW-bgQ5BptyB/view?usp=sharing	How do I use QuakePlate safely?	Who should NOT use QuakePlate?	What is QuakePlate and what does it do?					
How long is a Knee Pro session and what are the modes?	Charge device (first use 12+ hrs), mount device magnetically to wrap, place 3x5 thigh electrode and patella electrode, connect leads (red=patella). Start program and adjust intensity to a comfortable level. Quick Start: https://drive.google.com/file/d/17HlOn5mNzZxvJvtuPjW-bgQ5BptyB/view?usp=sharing	How do I set up Knee Pro?	How long is a Knee Pro session?	How often should Knee Pro be used?					
How often should Knee Pro be used?	Mode 1 (green) runs 30 min; Mode 2 (blue) runs 2 hours; unit shuts off at ~2.5 hours. Quick Start: https://drive.google.com/file/d/17HlOn5mNzZxvJvtuPjW-bgQ5BptyB/view?usp=sharing	What is NerveSpa Knee Pro?	How long is a Knee Pro session?	How often should Knee Pro be used?					
How do I care for Knee Pro electrodes and when to replace?	Recommended long-term program: 5 sessions/week, 130 sessions total (~325 hours over ~6 months); extend to 9–12 months for optimal results. Quick Start: https://drive.google.com/file/d/17HlOn5mNzZxvJvtuPjW-bgQ5BptyB/view?usp=sharing	How do I set up Knee Pro?	How often should Knee Pro be used?	What does Knee Pro feel like?					
What does Knee Pro therapy feel like?	Replace clear film after use; optionally unplug pigtails; each gel electrode typically lasts ~10–12 uses before replacement. Quick Start: https://drive.google.com/file/d/17HlOn5mNzZxvJvtuPjW-bgQ5BptyB/view?usp=sharing	What does Knee Pro feel like?	How long is a Knee Pro session and wh	How do I care for Knee Pro electrodes?					
Who should NOT use Knee Pro?	A mild, comfortable tingling. Reduce intensity if sharp/prickling, then increase to comfort.	How often should Knee Pro be used?	How do I set up Knee Pro?	What does Knee Pro feel like?					
What is NerveSpa Shoulder Pro?	Do not use with pacemakers, open wounds, when driving/operating equipment; consult physician for medical conditions. Quick Start: https://drive.google.com/file/d/17HlOn5mNzZxvJvtuPjW-bgQ5BptyB/view?usp=sharing	How long is a Knee Pro session and wh	How often should Knee Pro be used?	How do I care for Knee Pro electrodes and when to replace?					
How do I set up Shoulder Pro?	A shoulder wrap with Jstlm electrotherapy for regenerative therapy of the shoulder joint. Quick Start: https://drive.google.com/file/d/1QwjbWJhezrXTys7ha8bHr52qRlPjCf/view?usp=sharing	What is NerveSpa Knee Pro?	How do I set up Knee Pro?	What does Knee Pro feel like?					
How long is a Shoulder Pro session and what are the modes?	Charge device (first use 12+ hrs), place two 3" round electrodes on the wrap (leads to channel), connect red=top electrode, start program and adjust intensity. Quick Start: https://drive.google.com/file/d/1QwjbWJhezrXTys7ha8bHr52qRlPjCf/view?usp=sharing	How do I set up Shoulder Pro?	How long is a Shoulder Pro session?	How often should Shoulder Pro be used?					
How often should Shoulder Pro be used?	Mode 1 (green) runs 30 min; Mode 2 (blue) runs 2 hours; unit shuts off at ~2.5 hours. Quick Start: https://drive.google.com/file/d/1QwjbWJhezrXTys7ha8bHr52qRlPjCf/view?usp=sharing	How do I set up Shoulder Pro?	How long is a Shoulder Pro session and wh	How do I care for Shoulder Pro electrodes?					
How do I care for Shoulder Pro electrodes and when to replace?	Program: 5 sessions/week, 130 sessions (~325 hours over ~6 months); extend to 9–12 months depending on improvement. Quick Start: https://drive.google.com/file/d/1QwjbWJhezrXTys7ha8bHr52qRlPjCf/view?usp=sharing	How often should Shoulder Pro be used?	How do I set up Shoulder Pro?	What does Shoulder Pro feel like?					
What does Shoulder Pro therapy feel like?	Replace clear film after use; optionally unplug pigtails; each gel electrode typically lasts ~10–12 uses before replacement. Quick Start: https://drive.google.com/file/d/1QwjbWJhezrXTys7ha8bHr52qRlPjCf/view?usp=sharing	How long is a Shoulder Pro session and wh	How often should Shoulder Pro be u	How do I care for Shoulder Pro electrodes?					
Who should NOT use Shoulder Pro?	A mild, comfortable tingling. Reduce intensity if sharp/prickling, then increase to comfort.	How often should Shoulder Pro be u	How do I set up Shoulder Pro?	What does Shoulder Pro feel like?					
Is NerveSpa safe for seniors to use at home?	Contraindications include pacemakers/serious arrhythmia, pregnancy (especially abdominal), undiagnosed pain syndromes, epilepsy (consult doctor). Use only under physician supervision. Quick Start: https://drive.google.com/file/d/1QwjbWJhezrXTys7ha8bHr52qRlPjCf/view?usp=sharing	How long is a Shoulder Pro session	How often should Shoulder Pro be used?	How do I care for Shoulder Pro electrodes and when to replace?					
Can children use NerveSpa products?	Yes, NerveSpa systems are designed for home use under caregiver or self-monitoring. Seniors should ensure safe seating, proper water level (if footbath), and regular skin checks. Full User Guide: https://drive.google.com/file/d/1IEghxB0y0alwv3U4v9CDd7D4eulpxcgs/view?usp=sharing	Is NerveSpa safe?	What is NerveSpa Shoulder Pro?	How do I set up Shoulder Pro?					
How do I register my NerveSpa product for warranty?	NerveSpa products are not recommended for children under 18. Conductive garments and electrotherapy wraps specifically exclude pediatric use. Garments Guide: https://drive.google.com/file/d/18q1b0FWLzUbsCJmdrwuefh3oEnUSlC/view?usp=sharing	Who should NOT use NerveSpa Har	Can children use NerveSpa?	Is NerveSpa covered by insurance?					
What is included in the warranty?	Most devices include a manufacturer's warranty. Register your device online or via included card to activate. Keep proof of purchase. Shoulder Quick Start: https://drive.google.com/file/d/1QwjbWJhezrXTys7ha8bHr52qRlPjCf/view?usp=sharing	Who should NOT use the conductive	Is NerveSpa safe for seniors to use at h	What does treatment feel like?					
Who do I contact for technical support?	NerveSpa products include coverage against manufacturing defects for 1 year (standard). Consumables (salts, tablets, electrodes, sprays) are excluded. Shoulder Quick Start: https://drive.google.com/file/d/1QwjbWJhezrXTys7ha8bHr52qRlPjCf/view?usp=sharing	What is included in the warranty?	Who do I contact for technical support?	How long is the warranty period?					
What is the turnaround time for repairs or replacements?	Contact NerveSpa customer service at the official website or via phone/email provided in your manual. Manuals: https://drive.google.com/file/d/1IEghxB0y0alwv3U4v9CDd7D4eulpxcgs/view?usp=sharing	How do I register my NerveSpa prod	Who do I contact for technical support?	How long is the warranty period?					
What consumables come with NerveSpa systems?	Typical turnaround time is 2–3 weeks after receipt of your device at service center. Contact support for RMA. Manuals: https://drive.google.com/file/d/1IEghxB0y0alwv3U4v9CDd7D4eulpxcgs/view?usp=sharing	How do I register my NerveSpa prod	What is included in the warranty?	What is the turnaround time for repairs?					
How do I reorder consumables and accessories?	Hand & Footbath includes Epsom salts + efferesscent tablets; Conductive Garments require conductive spray; Knee/Shoulder include gel electrodes; QuakePlate requires no consumables. Quick Guide: https://drive.google.com/file/d/1ABWl8aU_hrnUSoRkmb3GoHl-rhbhTXlB/view?usp=sharing	Who do I contact for technical supp	What is included in the warranty?	How do I register my NerveSpa product for warranty?					
What accessories are included with each system?	Consumables (salts, tablets, conductive sprays, gel electrodes) and accessories (straps, cables, remotes) can be ordered through the NerveSpa website or authorized distributors. Manuals: https://drive.google.com/file/d/1IEghxB0y0alwv3U4v9CDd7D4eulpxcgs/view?usp=sharing	What consumables are used with the	How do I reorder consumables?	How often should consumables be replaced?					
What results can I expect from using NerveSpa Footbath?	Footbath includes leads/electrodes; Garments include gloves/socks; LED/Power Wraps include straps and power adapters; QuakePlate includes a remote; Knee/Shoulder include wraps and electrodes. Manuals: https://drive.google.com/file/d/1IEghxB0y0alwv3U4v9CDd7D4eulpxcgs/view?usp=sharing	What outcomes are expected from t	How often should the system be used?	What improvements do patients report?					
What results can I expect from using Conductive Garments?	Many users report improved circulation, reduced tingling, and greater comfort within weeks. Most consistent outcomes appear over 60–90 days of daily use. Full User Guide: https://drive.google.com/file/d/1IEghxB0y0alwv3U4v9CDd7D4eulpxcgs/view?usp=sharing	How often should conductive garmen	What does treatment feel like?	Who should NOT use the conductive garments?					
What results can I expect from using NerveBeam LED/Power Wraps?	Consistent daily sessions (4–6x/week) can reduce neuropathic pain and improve hand/foot comfort. Benefits often noted within 2–3 weeks, with stronger outcomes at 60–90 days. Garments Guide: https://drive.google.com/file/d/18q1b0FWLzUbsCJmdrwuefh3oEnUSlC/view?usp=sharing	How long is a NerveBeam session?	What are the benefits of red vs infrared?	Who should NOT use NerveBeam?					
What results can I expect from using QuakePlate?	Red light supports skin and surface tissue healing; red+infrared penetrates deeper for pain relief, inflammation reduction, and recovery. Outcomes may be seen after 2–4 weeks of consistent use. LED Wrap Manual: https://drive.google.com/file/d/11FLZ7E108Iyx4SdKSiBHP-NNskbOmC/view?usp=sharing	What is QuakePlate and what does i	How do I use QuakePlate safely?	Who should NOT use QuakePlate?					
What does treatment feel like?	Many users feel immediate muscle relaxation and reduced foot/leg tension. Circulation benefits may appear after 2–4 weeks of regular use. Manual: https://drive.google.com/file/d/1BNgb7AoSEeZTWMBtTh5R7r5xXs1b_n/view?usp=sharing	How often should Knee Pro be used?	How long is a Knee Pro session and wh	What does Knee Pro feel like?					
How quickly will I notice improvements?	Clinical protocol supports long-term improvement in OA/RA symptoms. Outcomes generally appear after 1–2 months of consistent use; full program (6–12 months) offers regenerative benefits. Knee Pro Quick Start: https://drive.google.com/file/d/17HlOn5mNzZxvJvtuPjW-bgQ5BptyB/view?usp=sharing	How often should Shoulder Pro be u	How long is a Shoulder Pro session and wh	What does Shoulder Pro therapy feel like?					
	Users often report reduced shoulder pain, increased mobility, and joint comfort after 1–2 months. Best outcomes occur with the full 6–12 month program. Shoulder Quick Start: https://drive.google.com/file/d/1QwjbWJhezrXTys7ha8bHr52qRlPjCf/view?usp=sharing	How long is a typical session?	What outcomes are expected from the H	Who should NOT use NerveSpa Hand & Footbath?					
	All NerveSpa Therapies should feel like a mild, comfortable tingling, warmth, or gentle vibration. Reduce intensity if prickling or burning occurs. Manuals: https://drive.google.com/file/d/1IEghxB0y0alwv3U4v9CDd7D4eulpxcgs/view?usp=sharing	What results can I expect from using	What outcomes are expected from the H	How often should the system be used?					
	Some users report relief within days; others require 2–6 weeks. Maximum benefits usually occur after 60–90 days of consistent daily use. Manuals: https://drive.google.com/file/d/1IEghxB0y0alwv3U4v9CDd7D4eulpxcgs/view?usp=sharing								

Question	Answer	Related Question 1	Related Question 2	Related Question 3		
Are results permanent?	Results depend on condition and compliance. Many users maintain improvements with continued regular use; discontinuing may allow symptoms to return gradually. Manuals: https://drive.google.com/file/d/1EgHxBOy0alwv3U4V9CDd7D4eulpxcgs/view?usp=sharing	What results can I expect from using How quickly will I notice improvements?	What outcomes are expected from the Hand & Foot Neuropathy System?			
Can NerveSpa be combined with other therapies?	Yes, NerveSpa is often used alongside medications, physical therapy, or lifestyle changes. Always consult your healthcare provider for combined treatment plans. Manuals: https://drive.google.com/file/d/1EgHxBOy0alwv3U4V9CDd7D4eulpxcgs/view?usp=sharing	Who should NOT use NerveSpa Har is NerveSpa safe for seniors to use at hc	What results can I expect from using Knee Pro?			
Can I use NerveSpa if I have a pacemaker?	Yes, NerveSpa is safe for use with pacemakers, defibrillators, or implanted stimulators without explicit physician clearance. Manuals: https://drive.google.com/file/d/1EgHxBOy0alwv3U4V9CDd7D4eulpxcgs/view?usp=sharing	Who should NOT use NerveSpa Har Who should NOT use Shoulder Pro?	Who should NOT use Knee Pro?			
Can I use NerveSpa during pregnancy?	NerveSpa systems are not recommended during pregnancy. This includes Footbath, Garments, Knee Pro, and Shoulder Pro. Manuals: https://drive.google.com/file/d/1EgHxBOy0alwv3U4V9CDd7D4eulpxcgs/view?usp=sharing	Who should NOT use NerveSpa Har Who should NOT use Shoulder Pro?	Is NerveSpa safe for seniors to use at home?			
Can NerveSpa be used if I have open wounds or skin conditions?	Do not use on areas with open wounds, rashes, infections, or broken skin. Wait until skin heals before resuming use. Manuals: https://drive.google.com/file/d/1EgHxBOy0alwv3U4V9CDd7D4eulpxcgs/view?usp=sharing	Who should NOT use NerveSpa Har Who should NOT use the conductive gar	Who should NOT use QuakePlate?			
Can I sleep with a NerveSpa device running?	No. NerveSpa devices (including NerveBeam wraps) should not be used while sleeping. Risk of burns, skin irritation, or unsafe postures may result. LED Wrap Manual: https://drive.google.com/file/d/11FLZ7E1O8iyxb4SdKSIBH-p-NNsKbOmC/view?usp=sharing	Who should NOT use NerveBeam?	How long is a NerveBeam session?	What does treatment feel like?		
Can NerveSpa products be used daily?	Yes. Most protocols recommend daily or near-daily use (4–6×/week) for best results. See product-specific guides. Footbath User Guide: https://drive.google.com/file/d/1EgHxBOy0alwv3U4V9CDd7D4eulpxcgs/view?usp=sharing	How often should the system be used	How often should Knee Pro be used?	How often should Shoulder Pro be used?		
Can NerveSpa help with arthritis?	Knee Pro and Shoulder Pro are specifically designed for osteoarthritis and rheumatoid arthritis support. Improvements in pain, mobility, and function are documented after consistent use. Knee Pro Quick Start: https://drive.google.com/file/d/17hLn0n5mNzxxvJvtUPJw-bgQ5BptyIB/view?usp=sharing	What is NerveSpa Knee Pro?	What is NerveSpa Shoulder Pro?	What results can I expect from using Knee Pro?		
Is NerveSpa HSA/FSA eligible?	Yes. NerveSpa devices can generally be purchased with HSA or FSA accounts. Check with your plan provider for documentation requirements. Manuals: https://drive.google.com/file/d/1EgHxBOy0alwv3U4V9CDd7D4eulpxcgs/view?usp=sharing	Is NerveSpa covered by insurance?	How do I register my NerveSpa product	What is included in the warranty?		
Do I need a prescription for NerveSpa?	Most products are available direct-to-consumer. Shoulder Pro may require a prescription in some jurisdictions as a regulated medical device. Shoulder Quick Start: https://drive.google.com/file/d/1QwkbWJhezzrXTys7ha8BhR5zqRJPJcF/view?usp=sharing	What is NerveSpa Shoulder Pro?	Is NerveSpa covered by insurance?	Who should NOT use Shoulder Pro?		
Can NerveSpa replace my medications?	No. NerveSpa is a complementary therapy. Continue medications unless advised otherwise by your physician. Manuals: https://drive.google.com/file/d/1EgHxBOy0alwv3U4V9CDd7D4eulpxcgs/view?usp=sharing	Can NerveSpa be combined with oth	What results can I expect from using Ner	What results can I expect from using Knee Pro?		
How do I clean and maintain the NerveSpa Footbath?	Empty basin after each session, dry electrodes and wires thoroughly, and store in a cool, dry place. Do not submerge controller. Quick Guide: https://drive.google.com/file/d/1ABWl8aU_hrvUSoRkmb3GoHl-rhbT7Xib/view?usp=sharing	What consumables are used with the Hand & Foot Neuropathy	Who should NOT use NerveSpa Hand & Footbath?			
How do I clean and maintain the conductive garments?	Hand-wash garments gently with mild soap, air-dry completely before reuse, and respray with conductive solution before each session. Garments Guide: https://drive.google.com/file/d/18q1bOFWVwLzUBscJmndruefhv3oEnUSlC/view?usp=sharing	What are NerveSpa Conductive Can	How often should conductive garments b	Who should NOT use the conductive garments?		
How do I clean and maintain NerveBeam wraps?	Unplug device, wipe gently with damp cloth, do not soak, and allow to fully dry before storage. Store flat and cool. LED Wrap Manual: https://drive.google.com/file/d/11FLZ7E1O8iyxb4SdKSIBH-p-NNsKbOmC/view?usp=sharing	How do I use the NerveBeam LED/P	Who should NOT use NerveBeam?	How long is a NerveBeam session?		
How do I store NerveSpa products between sessions?	Store devices in a cool, dry environment. Do not leave plugged in. Keep electrodes sealed on protective film. Manuals: https://drive.google.com/file/d/1EgHxBOy0alwv3U4V9CDd7D4eulpxcgs/view?usp=sharing	How do I clean and maintain the Ner	How do I care for Knee Pro electrodes a	How do I care for Shoulder Pro electrodes and when to replace?		
How often should consumables be replaced?	Exposm salt tablets: per session. Gel electrodes: ~10–12 uses. Conductive spray: reapply each session. Replace accessories when worn. Manuals: https://drive.google.com/file/d/1EgHxBOy0alwv3U4V9CDd7D4eulpxcgs/view?usp=sharing	What consumables come with Nerve	How do I reorder consumables and acce	What accessories are included with each system?		
How do I replace gel electrodes?	Peel gently from skin after each session, resal on plastic film, and replace after ~10–12 uses or if adhesion weakens. Shoulder Quick Start: https://drive.google.com/file/d/1QwkbWJhezzrXTys7ha8BhR5zqRJPJcF/view?usp=sharing	How do I care for Shoulder Pro elect	How do I care for Knee Pro electrodes a	How do I store NerveSpa products between sessions?		
How do I charge my NerveSpa device?	Plug supplied charger into wall outlet, connect to device port, and allow to fully charge (first use: 12+ hours). Recharge every 3–5 sessions depending on intensity. Manuals: https://drive.google.com/file/d/1EgHxBOy0alwv3U4V9CDd7D4eulpxcgs/view?usp=sharing	How should the battery be maintaine	How do I store NerveSpa products betw	What is included in the warranty?		
How do I troubleshoot if my device will not power on?	Ensure device is fully charged, hold power button for 2 seconds, check charger connections, and confirm outlet. If unresolved, contact support. Manuals: https://drive.google.com/file/d/1EgHxBOy0alwv3U4V9CDd7D4eulpxcgs/view?usp=sharing	How do I charge my NerveSpa devic	Who do I contact for technical support?	What is the turnaround time for repairs or replacements?		
How long is the warranty period?	Standard warranty is 1 year against manufacturing defects. Extended warranty options may be available through authorized distributors. Manuals: https://drive.google.com/file/d/1EgHxBOy0alwv3U4V9CDd7D4eulpxcgs/view?usp=sharing	What is included in the warranty?	How do I register my NerveSpa product	Who do I contact for technical support?		
Can I travel with my NerveSpa device?	Yes, devices are portable. Ensure charge is compatible with voltage at destination. Keep electrodes and consumables in sealed bags. Manuals: https://drive.google.com/file/d/1EgHxBOy0alwv3U4V9CDd7D4eulpxcgs/view?usp=sharing	How should the battery be maintaine	How do I store NerveSpa products betw	What accessories are included with each system?		
Are NerveSpa products FDA cleared?	NerveSpa devices are registered with the FDA as wellness/therapy devices. Check the NerveSpa website for current regulatory status by product. Manuals: https://drive.google.com/file/d/1EgHxBOy0alwv3U4V9CDd7D4eulpxcgs/view?usp=sharing	Is NerveSpa covered by insurance?	Is NerveSpa HSA/FSA eligible?	Do I need a prescription for NerveSpa?		
Do NerveSpa products come with instructions?	Yes, every product includes a printed quick start guide and an online manual. Manuals are also available here: https://drive.google.com/file/d/1EgHxBOy0alwv3U4V9CDd7D4eulpxcgs/view?usp=sharing	How do I register my NerveSpa prod	Who do I contact for technical support?	What is included in the warranty?		
How do I access NerveSpa demo videos?	Demo videos are available on the official YouTube channel: https://www.youtube.com/@nervespa and linked in the relevant FAQ answers.	What is NerveSpa (LED/Power Wrap	What are NerveSpa Conductive Garment	How do I use the NerveSpa Footbath?		
Can I share my NerveSpa device with others?	For hygiene reasons, consumables and electrodes should not be shared. If sharing, each user should use their own electrodes/sprays/salts. Manuals: https://drive.google.com/file/d/1EgHxBOy0alwv3U4V9CDd7D4eulpxcgs/view?usp=sharing	What consumables come with Nerve	How do I reorder consumables and acce	What accessories are included with each system?		
Are there clinical studies on NerveSpa?	NerveSpa technology is based on established electrotherapy and photobiomodulation principles. References are available on request via NerveSpa support. Manuals: https://drive.google.com/file/d/1EgHxBOy0alwv3U4V9CDd7D4eulpxcgs/view?usp=sharing	What results can I expect from using	What results can I expect from using NerveBeam LED/Power Wraps?			
What is the difference between LED Wrap and Power Wrap?	LED Wrap uses standard red/infrared output, while Power Wrap delivers higher energy for deeper tissue penetration and faster results. LED Wrap Manual: https://drive.google.com/file/d/11FLZ7E1O8iyxb4SdKSIBH-p-NNsKbOmC/view?usp=sharing Power Wrap Manual: https://drive.google.com/file/d/18q1bOFWVwLzUBscJmndruefhv3oEnUSlC/view?usp=sharing	What is NerveBeam (LED/Power Wr	How do I use the NerveBeam LED/Powe	What results can I expect from using NerveBeam LED/Power Wraps?		
What accessories require replacement most often?	Consumables like salts, effervescent tablets, conductive sprays, and gel electrodes need regular replacement. Straps and cables may also require periodic replacement. Manuals: https://drive.google.com/file/d/1EgHxBOy0alwv3U4V9CDd7D4eulpxcgs/view?usp=sharing	How often should consumables be r	How do I reorder consumables and acce	What accessories are included with each system?		
What is included in the starter kit for each product?	Each system includes device, leads/cables, and starter consumables (varies by product). Manuals list exact contents. Manuals: https://drive.google.com/file/d/1EgHxBOy0alwv3U4V9CDd7D4eulpxcgs/view?usp=sharing	What accessories are included with	What consumables come with NerveSpa	How do I reorder consumables and accessories?		
How do I know which NerveSpa product is right for me?	Footbath is best for peripheral neuropathy in feet; Garments for hands/feet neuropathy; LED/Power Wraps for localized pain/inflammation; Knee Pro/Shoulder Pro for arthritis; QuakePlate for vibration therapy and circulation. Manuals: https://drive.google.com/file/d/1EgHxBOy0alwv3U4V9CDd7D4eulpxcgs/view?usp=sharing	What results can I expect from using	What is NerveSpa Knee Pro?	What is QuakePlate and what does it do?		
Can I combine multiple NerveSpa therapies?	Yes, many users combine devices (e.g., Footbath + LED Wrap, or Knee Pro + QuakePlate) for complementary effects. Consult physician before combining electrotherapy and phototherapy. Manuals: https://drive.google.com/file/d/1EgHxBOy0alwv3U4V9CDd7D4eulpxcgs/view?usp=sharing	Can NerveSpa replace my medicati	Can NerveSpa be combined with other t	What results can I expect from using Knee Pro?		
What are the benefits of red vs infrared light?	Red (visible) light promotes skin healing; infrared penetrates deeper to aid joints, muscles, and circulation. Many wraps combine both for maximum effect. LED Wrap Manual: https://drive.google.com/file/d/11FLZ7E1O8iyxb4SdKSIBH-p-NNsKbOmC/view?usp=sharing	What results can I expect from using	How long is a NerveBeam session?	Who should NOT use NerveBeam?		
What is the difference between NerveSpa Footbath and Conductive Garments?	Footbath provides immersion therapy with salts and electrodes for feet; Garments use gloves/socks with conductive spray for flexible use (hands/feet). Garments Guide: https://drive.google.com/file/d/18q1bOFWVwLzUBscJmndruefhv3oEnUSlC/view?usp=sharing	What are NerveSpa Conductive Can	How do I use the NerveSpa Footbath?	What results can I expect from using NerveSpa Footbath?		
Can I use QuakePlate together with Footbath?	No. Use QuakePlate separately from Footbath sessions. Combining could cause unsafe vibration with water/electrodes. QuakePlate Manual: https://drive.google.com/file/d/1BNgb7AoSEeZTWmbTTh5r7rx5x1b_n/view?usp=sharing	What is QuakePlate and what does	How do I use QuakePlate safely?	Who should NOT use QuakePlate?		
Can I use QuakePlate before or after Knee/Shoulder Pro sessions?	Yes, QuakePlate can be used before Knee/Shoulder therapy to relax muscles, or after sessions to enhance circulation. Allow a break between devices. QuakePlate Manual: https://drive.google.com/file/d/1BNgb7AoSEeZTWmbTTh5r7rx5x1b_n/view?usp=sharing	What results can I expect from using	How often should Knee Pro be used?	How often should Shoulder Pro be used?		
How do I know if my electrodes are still effective?	Good electrodes should stick well and conduct current evenly. Replace when tingling is uneven, adhesion weakens, or after ~10–12 uses. Shoulder Quick Start: https://drive.google.com/file/d/1QwkbWJhezzrXTys7ha8BhR5zqRJPJcF/view?usp=sharing	How do I care for Shoulder Pro elect	How do I care for Knee Pro electrodes a	What if stimulation feels weak?		
Can I use NerveSpa with compression socks or braces?	No. Electrodes should contact bare skin. Do not layer NerveSpa therapy with compression garments or braces unless instructed by a provider. Manuals: https://drive.google.com/file/d/1EgHxBOy0alwv3U4V9CDd7D4eulpxcgs/view?usp=sharing	Who should NOT use NerveSpa Har	Who should NOT use the conductive gar	What does treatment feel like?		
Are there side effects from NerveSpa?	Most users tolerate therapy well. Mild redness, tingling, or skin irritation may occur; usually resolves quickly. Discontinue if persistent discomfort. Manuals: https://drive.google.com/file/d/1EgHxBOy0alwv3U4V9CDd7D4eulpxcgs/view?usp=sharing	Who should NOT use Knee Pro?	Who should NOT use Shoulder Pro?	Who should NOT use NerveBeam?		
What precautions should I take before starting therapy?	Read the Quick Start guide, ensure skin is clean/dry, use recommended consumables, and start with lower intensity. Increase gradually as tolerated. Quick Guide: https://drive.google.com/file/d/1ABWl8aU_hrvUSoRkmb3GoHl-rhbT7Xib/view?usp=sharing	What does treatment feel like?	Who should NOT use NerveSpa Hand &	What improvements do patients report?		
Can NerveSpa devices be used in clinics or only at home?	NerveSpa systems are designed for home use but are also used in physical therapy and pain clinics. Check with your provider for supervised protocols. Manuals: https://drive.google.com/file/d/1EgHxBOy0alwv3U4V9CDd7D4eulpxcgs/view?usp=sharing	Is NerveSpa safe for seniors to use	Do I need a prescription for NerveSpa?	Are there clinical studies on NerveSpa?		
Can I return a NerveSpa product if it doesn't work for me?	Return policies depend on distributor or seller. Typically 30-day return windows apply if device is undamaged. Check seller terms. Manuals: https://drive.google.com/file/d/1EgHxBOy0alwv3U4V9CDd7D4eulpxcgs/view?usp=sharing	What is included in the warranty?	Who do I contact for technical support?	How do I register my NerveSpa product for warranty?		
What conditions does NerveSpa support?	NerveSpa is designed for neuropathy, arthritis, chronic pain, and circulation issues. Results vary by condition and adherence. Manuals: https://drive.google.com/file/d/1EgHxBOy0alwv3U4V9CDd7D4eulpxcgs/view?usp=sharing	What results can I expect from using	Can NerveSpa help with arthritis?	What improvements do patients report?		
Are NerveSpa products safe to use with medications?	Yes, NerveSpa can be used while continuing most medications. Always consult a physician before altering prescriptions. Manuals: https://drive.google.com/file/d/1EgHxBOy0alwv3U4V9CDd7D4eulpxcgs/view?usp=sharing	Can NerveSpa replace my medicati	Can NerveSpa be combined with other t	What precautions should I take before starting therapy?		
What should I do if I feel discomfort during a session?	Stop immediately, reduce intensity, and recheck electrode placement or consumables. Resume only if comfortable. Consult a physician if persistent. Manuals: https://drive.google.com/file/d/1EgHxBOy0alwv3U4V9CDd7D4eulpxcgs/view?usp=sharing	How do treatment feel like?	What if stimulation feels weak?	Are there side effects from NerveSpa?		
Do NerveSpa devices need calibration?	No, all devices are pre-calibrated. Only consumables/electrodes require replacement as instructed. Manuals: https://drive.google.com/file/d/1EgHxBOy0alwv3U4V9CDd7D4eulpxcgs/view?usp=sharing	How do I replace gel electrodes?	How do I charge my NerveSpa device?	How do I troubleshoot if my device will not power on?		
How long do NerveSpa devices typically last?	With proper care, most units last 3–5 years. Consumables require periodic replacement. Manuals: https://drive.google.com/file/d/1EgHxBOy0alwv3U4V9CDd7D4eulpxcgs/view?usp=sharing	How do I clean and maintain the Ner	How do I store NerveSpa products betw	What accessories require replacement most often?		
Where can I find NerveSpa product manuals?	All current manuals are available in print with each product and online here: https://drive.google.com/file/d/1EgHxBOy0alwv3U4V9CDd7D4eulpxcgs/view?usp=sharing	Do NerveSpa products come with in	How do I access NerveSpa demo videos	Who do I contact for technical support?		
Where can I watch NerveSpa YouTube videos?	Visit the official YouTube channel: https://www.youtube.com/@nervespa for demos, setup guidance, and educational content.	How do I access NerveSpa demo vi	What results can I expect from using	How do I use the NerveBeam LED/Power Wrap?		
Can I get medical advice from NerveSpa?	No. NerveSpa provides device guidance only. For diagnosis, prescriptions, or clinical decisions, consult your healthcare provider. Manuals: https://drive.google.com/file/d/1EgHxBOy0alwv3U4V9CDd7D4eulpxcgs/view?usp=sharing	Are NerveSpa products safe to use	Can NerveSpa replace my medications?	What precautions should I take before starting therapy?		
Show me a demo	Which product would you like a demo for? Options: Footbath • Conductive Garments • LED Wrap • Power Wrap • QuakePlate • Knee Pro • Shoulder Pro. You can also say "all" to see demos/manuals for every device. Footbath — Demos: https://www.youtube.com/@nervespa • Manuals: Quick Guide https://drive.google.com/file/d/1ABWl8aU_hrvUSoRkmb3GoHl-rhbT7Xib/view?usp=sharing User Guide https://drive.google.com/file/d/1EgHxBOy0alwv3U4V9CDd7D4eulpxcgs/view?usp=sharing Conductive Garments — Demos: https://www.youtube.com/@nervespa • Manual: https://drive.google.com/file/d/18q1bOFWVwLzUBscJmndruefhv3oEnUSlC/view?usp=sharing LED Wrap — Demos: https://www.youtube.com/@nervespa • Manual: https://drive.google.com/file/d/11FLZ7E1O8iyxb4SdKSIBH-p-NNsKbOmC/view?usp=sharing Power Wrap — Demos: https://www.youtube.com/@nervespa • Manual: https://drive.google.com/file/d/18q1bOFWVwLzUBscJmndruefhv3oEnUSlC/view?usp=sharing QuakePlate — Demos: https://www.youtube.com/@nervespa • Manual: https://drive.google.com/file/d/1BNgb7AoSEeZTWmbTTh5r7rx5x1b_n/view?usp=sharing Knee Pro — Demos: https://www.youtube.com/@nervespa • Quick Start: https://drive.google.com/file/d/17hLn0n5mNzxxvJvtUPJw-bgQ5BptyIB/view?usp=sharing Shoulder Pro — Demos: https://www.youtube.com/@nervespa • Quick Start: https://drive.google.com/file/d/1QwkbWJhezzrXTys7ha8BhR5zqRJPJcF/view?usp=sharing	Show me demos for all devices	Show me a demo for Footbath	Show me a demo for Knee Pro		
Show me demos for all devices	YouTube demos: https://www.youtube.com/@nervespa • Manuals: Quick Guide https://drive.google.com/file/d/1ABWl8aU_hrvUSoRkmb3GoHl-rhbT7Xib/view?usp=sharing User Guide https://drive.google.com/file/d/1EgHxBOy0alwv3U4V9CDd7D4eulpxcgs/view?usp=sharing	Show me a demo	Show me demos for all devices	How should the Hand & Foot Neuropathy System be set up?		
Show me a demo for Footbath	YouTube demos: https://www.youtube.com/@nervespa • Manual: https://drive.google.com/file/d/18q1bOFWVwLzUBscJmndruefhv3oEnUSlC/view?usp=sharing	Show me a demo	Show me demos for all devices	How do I use the conductive gloves or socks?		
Show me a demo for Conductive Garments	YouTube demos: https://www.youtube.com/@nervespa • Manual: https://drive.google.com/file/d/11FLZ7E1O8iyxb4SdKSIBH-p-NNsKbOmC/view?usp=sharing	Show me a demo	Show me demos for all devices	How do I use the NerveBeam LED/Power Wrap?		
Show me a demo for LED Wrap	YouTube demos: https://www.youtube.com/@nervespa • Manual: https://drive.google.com/file/d/11FLZ7E1O8iyxb4SdKSIBH-p-NNsKbOmC/view?usp=sharing	Show me a demo	Show me demos for all devices	How do I use the NerveBeam LED/Power Wrap?		
Show me a demo for Power Wrap	YouTube demos: https://www.youtube.com/@nervespa • Manual: https://drive.google.com/file/d/1WdJ7d9KnGU25nkKp-T7o9aupF5/view?usp=sharing	Show me a demo	Show me demos for all devices	What is QuakePlate and what does it do?		
Show me a demo for QuakePlate	YouTube demos: https://www.youtube.com/@nervespa • Manual: https://drive.google.com/file/d/1BNgb7AoSEeZTWmbTTh5r7rx5x1b_n/view?usp=sharing	Show me a demo	Show me demos for all devices	How do I set up Knee Pro?		
Show me a demo for Knee Pro	YouTube demos: https://www.youtube.com/@nervespa • Quick Start: https://drive.google.com/file/d/17hLn0n5mNzxxvJvtUPJw-bgQ5BptyIB/view?usp=sharing	Show me a demo	Show me demos for all devices	How do I set up Shoulder Pro?		
Show me a demo for Shoulder Pro	YouTube demos: https://www.youtube.com/@nervespa • Quick Start: https://drive.google.com/file/d/1QwkbWJhezzrXTys7ha8BhR5zqRJPJcF/view?usp=sharing					

Question	Answer	Related Question 1	Related Question 2	Related Question 3				
Troubleshoot Knee Pro	Charge fully (first use 12+ hrs). Hold power 2s. Check lead connections (red--patella), confirm Mode (green=30 min, blue=2 hrs). Ensure gel electrodes adhere; replace after ~10–12 uses. Increase intensity gradually; hydrate skin. Auto shut-off at ~2.5 hrs is normal. Quick Start: https://drive.google.com/file/d/17hL0n5mNz2xvJvUjW-bq55p8tB/view?usp=sharing	How do I set up Knee Pro?	How long is a Knee Pro session and what are the modes?	How do I care for Knee Pro electrodes and when to replace?				
Troubleshoot Shoulder Pro	Charge fully (first use 12+ hrs). Hold power 2s. Check leads (red--top electrode). Verify Mode (green=30 min, blue=2 hrs). Ensure gel electrodes adhere; replace after ~10–12 uses. Adjust intensity to mild, comfortable tingle. Auto shut-off at ~2.5 hrs is normal. Quick Start: https://drive.google.com/file/d/1QwkbWJheezrXTys7ha8bHr52qRlPjCf/view?usp=sharing	How do I set up Shoulder Pro?	How long is a Shoulder Pro session and what are the modes?	How do I care for Shoulder Pro electrodes and when to replace?				
How do I clean and store Knee Pro?	Wipe wrap with a lightly damp cloth; do not soak. Keep ports dry. Reseal gel electrodes on protective film; replace after ~10–12 uses. Store cool/dry; unplug pigtails if desired. Quick Start: https://drive.google.com/file/d/17hL0n5mNz2xvJvUjW-bq55p8tB/view?usp=sharing	How do I care for Knee Pro electrodes and when to replace?	Troubleshoot Knee Pro	How often should Knee Pro be used?				
How do I clean and store Shoulder Pro?	Wipe wrap with a lightly damp cloth; do not soak. Keep ports dry. Reseal gel electrodes on protective film; replace after ~10–12 uses. Store cool/dry; unplug pigtails if desired. Quick Start: https://drive.google.com/file/d/1QwkbWJheezrXTys7ha8bHr52qRlPjCf/view?usp=sharing	How do I care for Shoulder Pro electrodes and when to replace?	Troubleshoot Shoulder Pro	How often should Shoulder Pro be used?				
How long is a QuakePlate session?	Start with 10 minutes and increase up to ~20 minutes as tolerated. Monitor comfort; stop if discomfort occurs. Manual: https://drive.google.com/file/d/1BNgb7JAoSEezTWMbTTh5R7r5xS1b_n/view?usp=sharing	What is QuakePlate and what does it do?	How do I use QuakePlate safely?	Who should NOT use QuakePlate?				
How often should QuakePlate be used?	Use 1–2x daily or per clinic protocol on non-irritated skin. Do not combine with water/electrotherapy sessions. Manual: https://drive.google.com/file/d/1BNgb7JAoSEezTWMbTTh5R7r5xS1b_n/view?usp=sharing	What is QuakePlate and what does it do?	What is QuakePlate and what does it do?	How do I clean and maintain QuakePlate?				
Welcome (Clinic)	Welcome to NerveSpa Clinic. Choose a section or type a question.	Setup	Protocols	Troubleshoot				
Main Menu (Clinic)	Setup Protocols Troubleshoot Billing Orders Warranty	Billing	Orders	Warranty				
Choose a device	Select a device to continue: Footbath Conductive Garments LED Wrap Power Wrap QuakePlate Knee Pro Shoulder Pro	Protocols	Orders	Troubleshoot				
How do I clean and maintain the conductive garments?	Hand-wash gloves/socks gently with mild soap; do not wring; air-dry completely before reuse. Spray inside/outside with conductive solution before each session. Store cool/dry. Manual: https://drive.google.com/file/d/18q1b0FWvLUZuBscJmdwuefh30cEnUSlcv/view?usp=sharing	What are NerveSpa Conductive Garments (gloves/socks)?	How often should conductive garments be used?	How do I use the conductive gloves or socks?				
How often should the LED Wrap be used?	Use once or twice daily on intact skin as tolerated; stop if irritation occurs. One session is ~20 minutes (auto-shutoff). Follow clinic protocol. Manual: https://drive.google.com/file/d/11FtLZ7E1O8lyxb4dSkSIBHn-NNskbOmCv/view?usp=sharing	How long is a NerveBeam session?	How do I use the NerveBeam LED/Power Wrap?	Who should NOT use NerveBeam?				
How often should the Power Wrap be used?	Use once or twice daily on intact skin as tolerated; stop if irritation occurs. One session is ~20 minutes (auto-shutoff). Follow clinic protocol. Manual: https://drive.google.com/file/d/1wWDJ7dD9kanGU925nk0Kp-T7o9faupF5/view?usp=sharing	How long is a NerveBeam session?	How do I use the NerveBeam LED/Power Wrap?	Who should NOT use NerveBeam?				
Troubleshoot QuakePlate	Check power is connected to 110V AC, switch is ON, and connections are dry. If remote is unresponsive, replace batteries and re-pair per manual. Stop use with any adverse reaction. Manual: https://drive.google.com/file/d/1BNgb7JAoSEezTWMbTTh5R7r5xS1b_n/view?usp=sharing	What is QuakePlate and what does it do?	How do I use QuakePlate safely?	How do I clean and maintain QuakePlate?				
Show videos	Which product would you like videos for? Options: Footbath • Conductive Garments • LED Wrap • Power Wrap • QuakePlate • Knee Pro • Shoulder Pro. You can also say "all" for every device.	Show me a demo	Show me demos for all devices	Where can I watch NerveSpa YouTube videos?				
Product video	Sure—pick a product, or say "all" to see every device: Footbath • Conductive Garments • LED Wrap • Power Wrap • QuakePlate • Knee Pro • Shoulder Pro. Providers can reach out for support in the following ways:	Show me a demo	Show me demos for all devices	Where can I watch NerveSpa YouTube videos?				
contact info	• Chatbot – Use the "Contact Support" option directly in the chatbot.							
phone number	• Email – support@nervespa.com							
email id	• Phone – 1-800-239-7880 (available during business hours)							
	• Provider Portal – coming soon							
	Phone – 1-800-239-7880 (available during business hours)							
	Email – support@nervespa.com							